



Analysis of Cyromazine in Poultry Feed Using a QuEChERS Approach

UCT Products:

ECMSSA50CT-MP (6 g anhydrous MgSO₄ and 1.5 g Na Acetate)

EEC18156 (500 mg endcapped C18, 6 mL cartridge)

October 12, 2010

Introduction

This summary outlines a QuEChERS procedure for the analysis of the insecticide cyromazine (Trigard or Larvadex) in poultry feed by LC-MS/MS. Processing time is significantly faster than EPA method AG-555 and uses less solvent. Modifications include adding glacial acetic acid to the acetonitrile to increase extraction efficiency.

Procedure

1) Sample Preparation

- a) Homogenize 2 grams of poultry feed and add to a 50 mL centrifuge tube
- b) Add 10 mLs of acetonitrile/acetic acid (75:25)
- c) Sonicate at 50/60 Hz for 15 minutes
- d) Add the contents of **ECMSSA50CT-MP** pouch and shake for 1 minute
- e) Centrifuge at 3400 rpm for 10 minutes
- f) Transfer 1 mL of supernatant to a calibrated test tube and add 9 mL of water: acetonitrile (95:5) with 0.1% acetic acid

2) Sample Cleanup

- a) Add the 10 mLs from 1) f) above to a **EEC18156** cartridge and elute dropwise
- b) Filter eluant using a 0.45 µm Teflon filter (Millipore, Billerica, MA) or equivalent
- c) Transfer 2 mL of eluant to an HPLC vial for analysis by LC-MS/MS

3) Analysis LC-MS/MS

- Waters Alliance 2695 HPLC (Waters) coupled with a micromass Quattro Micro triple-quadrupole mass spectrometer (Micromass, Manchester, U.K.) or equivalent
- **HPLC conditions:**
 - **Guard column** (Alltima, C18, 5 μ m, 2.1 x 7.5 mm, Deerfield, IL) or equivalent
 - **Analytical column** (Alltima, C18, 5 μ m, 2.1 x 250 mm, Waters) or equivalent
 - **Mobile phase:** (A) acetonitrile with 0.1% formic acid and (B) water with 0.1% formic acid
 - **Gradient:**
 - 0-2 min, 5%A
 - 2-5 min from 5 to 10% A
 - 5-5.5 min from 10 to 90% A
 - 5.5-8 min 90 to 5% A
 - 8-10 min, from 90 to 5% A
 - 10-12 min, 5% A
 - **Flow rate** 0.2 mL/minute
 - **Injection volume:** 25 μ L
- **Mass Spectrometer**
 - Positive ion mode electrospray ionization
 - Monitor the ion transition of the parent ion (m/z 167) to the product ion (m/z 85) in multiple reaction monitoring (MRM)

Mass Spectrometry Conditions for Cyromazine Quantitation

capillary voltage	3.1 kV
cone voltage	65 V
collision energy	21-24 V
source temperature	120° C
desolvation temperature	350° C
cone gas flow	135 L/h
desolvation gas flow rate	750 L/h
collision gas	Argon
parent ion	(m/z) 167
product ion	(m/z) 85

*Summarized with permission from Xia, Kang, Atkins, Jack et al, "Analysis of Cyromazine in Poultry Feed Using the QuEChERS Method Coupled with LC-MS/MS" J. Agric. Food Chem, DOI:10.1021/jf9034282

Listing of instrument manufacturers does not constitute endorsement by UCT

DCN-012101-196