

# Antidepressants & Painkillers In Blood and Urine Using: 200 mg Clean Screen® Extraction Column



## UCT Part Numbers

**ZSDAU020**  
Clean Screen® DAU  
10 mL, 200 mg sorbent  
Without Tips

## Procedure:

### 1. Prepare Sample

- To 1 mL of 100 mM phosphate buffer (pH= 6) add internal standard.\*
- Add 1 mL of blood or urine. Add 2 mL of 100 mM phosphate buffer (pH= 6). Mix/vortex.
- Sample pH should be  $6.0 \pm 0.5$ .
- Adjust pH accordingly with 100 mM monobasic or dibasic sodium phosphate.
- Mix/vortex.
- Centrifuge as appropriate.

### 2. Condition Clean Screen® Extraction Column

- 1 x 3 mL MeOH.
  - 1 x 3 mL H<sub>2</sub>O.
  - 1 x 1 mL 100 mM phosphate buffer (pH= 6.0).
- Note:** Aspirate at < 3 inches Hg to prevent sorbent drying.

### 3. Apply Sample

- Load at 1-2 mL/minute.

### 4. Wash Column

- 1 x 3 mL D.I. H<sub>2</sub>O.
- 1 x 3 mL 1% acetic acid.
- 1 x 3 mL Methanol.
- Dry column (5 minutes at >10 inches Hg).

### 5. Elute Antidepressants/Painkillers

- 1 x 3 mL ethyl acetate; acetonitrile: ammonia (78: 20: 2 v/v)
- Collect eluate at 1-2 mL /minute.

### 6. Evaporation

- Evaporate eluate under a gentle stream of nitrogen < 40 °C.

### 7. Reconstitute

- Reconstitute sample in 100 µL of Methanol.
- Inject 5 µL.

## Instrument Conditions

|                    |                                                    |
|--------------------|----------------------------------------------------|
| Column             | 150 x 2.1 mm (3.5 µm) Zorbax: Agilent Technologies |
| Mobile Phase       | Acetonitrile: 0.1% Formic acid (33:67)             |
| Column Temperature | ambient                                            |
| Flow Rate          | 0.35 mL/minute                                     |
| Detector           | API 2000 MS/MS                                     |



| Compound            | MRM Transition | Cerilliant #  |
|---------------------|----------------|---------------|
| Amitriptyline       | 278.8/91.1     | A-923         |
| *Amitriptyline-D3   | 281.2/91.2     | B-130284 N-10 |
| Diphenhydramine     | 256.2/167.1    | D-015         |
| *Diphenhydramine-D3 | 259.2/167.1    | D-017         |
| Doxepin             | 280.2/107.1    | D-005         |
| EDDP                | 278.2/234.2    | E-022         |
| *EDDP-D3            | 281.4/234.3    | E-021         |
| Methadone           | 310.2/105.1    | M-007         |
| *Methadone-D9       | 319.2/268.3    | M-088         |
| Nortriptyline       | 264.2/91.1     | N-907         |
| Norpropoxyphene     | 326.2/44.1     | N-913         |
| *Norpropoxyphene-D5 | 331.1/267.1    | N-904         |
| Propoxyphene        | 340.2/58.1     | P-011         |
| *Propoxyphene-D11   | 351.3/64.0     | P-013         |
| Sertraline          | 308.1/161.0    | S-006         |
| Tramadol            | 264.2/58.1     | T-027         |
| *Tramadol-D3        | 268.2/58.0     | T-029         |
| Venlafaxine         | 278.2/58.2     | V-004         |
| Zolpidem            | 308.2/235      | Z-901         |

UCT, LLC • 2731 Bartram Road • Bristol, PA 19007 800.385.3153 • 215.781.9255

[www.unitedchem.com](http://www.unitedchem.com) Email: [methods@unitedchem.com](mailto:methods@unitedchem.com)

©UCT, LLC 2008 • All rights reserved

