

Fentanyl-Norfentanyl On Oral Swabs Using: 200 mg Clean Screen® Extraction Column



UCT Part Numbers

ZSDAU020
Clean Screen® DAU
10 mL, 200 mg sorbent
Without Tips

Procedure:

1. Prepare Sample

Preparation of Standards:

- To separate tube add 0, 1, 5, 10 ng of Fentanyl / Norfentanyl in methanol.
- Evaporate off the solvent. Add 100 µL of drug free oral fluid. Vortex mix and allow to stand for 30 minutes.
- Take clean, dry (drug free) swab and swab up the oral fluid and allow standing for 15 minutes. Remove oral swab.

Sample Pretreatment:

- To 200 µL of methanol (pH 6) add internal standard.*
- Insert oral swab into methanol and mix for 1 minute, add a further 100 µL of methanol, allow to stand for 10 minutes.
- Remove swab and 3 mL of 100 mM phosphate buffer (pH 6). Vortex and centrifuge as appropriate.

2. Condition Column

- 1 x 3 mL MeOH.
- 1 x 3 mL H₂O.
- 1 x 1 mL 100 mM phosphate buffer (pH= 6).

Note: Aspirate at < 3 inches Hg to prevent sorbent drying.

3. Apply Sample

- Load at 1-2 mL/minute.

4. Wash Column

- 1 x 3 mL D.I. H₂O.
- 1 x 3 mL 1% actate acid.
- 1 x 3 mL Methanol.
- Dry column (5 minutes at > 10 inches Hg).

5. Elute Fentanyl/Norfentanyl

- 1 x 3 mL ethyl acetate; acetonitrile: ammonia (78: 20: 2 v/v)
- Collect eluate at 1-2 mL /minute.

6. Evaporation

- Evaporate to dryness at < 40 °C.

7. Reconstitute sample in 20 µL of methanol.


- Inject 5 µL.



Instrument Conditions

Column	150 x 2.1 mm (3 µm) PFP Thermofisher
Mobile Phase	Acetonitrile: 0.1% Formic acid (60:40)
Flowrate	0.35 mL/minute
Column Temperature	ambient
Detector	API 2000 MS/MS

Compound	MRM Transition	Cerilliant #
Fentanyl	333.2/188.3	F-002
*Fentanyl-D5	342.2/188.2	F-001
Norfentanyl	233.2/84.1	N-031
*Norfentanyl-D5	238.3/84.1	N-030



Recovery: > 90% (N=100)

LOD: 1 ng/m

UCT, LLC • 2731 Bartram Road • Bristol, PA 19007 800.385.3153 • 215.781.9255

www.unitedchem.com Email: methods@unitedchem.com

©UCT, LLC 2008 • All rights reserved

