

Lysergic Acid Diethylamide (LSD) In Urine For GC or GC/MS Confirmations Using: 200 mg Clean Screen® Extraction Column



UCT Part Numbers

ZSDAU020 Clean Screen® DAU 10 mL, 200 mg sorbent Without Tips	Or	ZCDAU020 Clean Screen® DAU 10 mL, 200 mg sorbent With CLEAN-THRU® Tips
---	----	--

Procedure:

1. Prepare Sample

- To 4 mL of D.I. H₂O internal standard* add 1 mL of blood or plasma/serum
- Mix/vortex and let stand 5 minutes.
- Centrifuge for 10 minutes at 2000 rpm and discard pellet.
- Add 2 mL 100 mM phosphate buffer (pH = 6.0).
- Mix/vortex. Sample pH should be 6.0 ± 0.5.
- Adjust pH accordingly with 100 mM monobasic or dibasic sodium phosphate.

2. Condition Clean Screen® Extraction Column

- 1 x 3 mL CH₃OH.
- 1 x 3 mL D.I. H₂O.
- 1 x 1 mL 100 mM phosphate buffer (pH = 6.0).

Note: Aspirate at < 3 inches Hg. to prevent sorbent drying.

3. Apply Sample

- Load at 1 mL/minute.

4. Wash Column

- 1 x 3 mL D.I. H₂O.
- 1 x 1 mL 100 mM acetic acid.
- 1 x 3 mL CH₃OH.
- Dry column (5 minutes at > 10 inches Hg).

5. Elute LSD

- 1 x 3 mL CH₂Cl₂/IPA/NH₄OH (78:20:2);
 - Collect eluate at 1 mL/minute.
- Note:** Prepare elution solvent daily. Add IPA/NH₄OH, mix, then add CH₂Cl₂ (pH 11-12).

6. Dry Eluate

- Evaporate to dryness at < 40°C.

7. Derivatize

- Add 20 µL ethyl acetate and 20 µL BSTFA (with 1% TMCS)***.
 - Overlay with N₂ and cap. Mix/vortex.
 - React 30 minutes at 70°C.
- Remove from heat source to cool.

Note: Do not evaporate BSTFA solution.

8. Quantitate

- Inject 1 to 2 µL onto gas chromatograph.
- For MSD monitor the following ions:

Analyte	Primary Ion***	Secondary	Tertiary	Cerilliant #
LSD-D ₃ -TMS*	298	296	271	L-006
LSD-TMS	395	293	268	L-005

* Suggested internal standard for GC/MS: D₃-LSD

*** Part # SBSTFA-1-1,10,25,100

**** Quantitation ion



UCT, LLC • 2731 Bartram Road • Bristol, PA 19007 800.385.3153 • 215.781.9255
www.unitedchem.com Email: methods@unitedchem.com
©UCT, LLC 2008 • All rights reserved

