


# Analysis of Cyromazine in Poultry Feed Using a QuEChERS Approach



## UCT Part Numbers

**ECMSSA50CT-MP**  
6 g anhydrous MgSO<sub>4</sub>  
and 1.5 g Na Acetate

**EEC18156**  
500 mg endcapped C18,  
6 mL cartridge

## Introduction:

This summary outlines a QuEChERS procedure for the analysis of the insecticide cyromazine (Trigard or Larvadex) in poultry feed by LC-MS/MS. Processing time is significantly faster than EPA method AG-555 and uses less solvent. Modifications include adding glacial acetic acid to the acetonitrile to increase extraction efficiency.

## Procedure:

### 1. Sample Preparation

- Homogenize 2 grams of poultry feed and add to a 50 mL centrifuge tube
- Add 10 mLs of acetonitrile/acetic acid (75:25)
- Sonicate at 50/60 Hz for 15 minutes
- Add the contents of **ECMSSA50CT-MP** pouch and shake for 1 minute
- Centrifuge at 3400 rpm for 10 minutes
- Transfer 1 mL of supernatant to a calibrated test tube and add 9 mL of water: acetonitrile (95:5) with 0.1% acetic acid

### 2. Sample Cleanup


- Add the 10 mLs from 1) f) above to a **EEC18156** cartridge and elute dropwise
- Filter eluant using a 0.45 µm Teflon filter (Millipore, Billerica, MA) or equivalent
- Transfer 2 mL of eluant to an HPLC vial for analysis by LC-MS/MS

### 3. Analysis LC-MS/MS

- Waters Alliance 2695 HPLC (Waters) coupled with a micromass Quattro Micro triplequadrupole mass spectrometer (Micromass, Manchester, U.K.) or equivalent

#### HPLC Conditions

<b>Guard column</b>	(Alltima, C18, 5 µm, 2.1 x 7.5 mm, Deerfield, IL) or equivalent
<b>Analytical column</b>	(Alltima, C18, 5 µm, 2.1 x 250 mm, Waters) or equivalent
<b>Flow rate</b>	0.2 mL/minute
<b>Injection volume</b>	25 µL




Gradient Program		
Time (min)	(A) acetonitrile with 0.1% formic acid	(B) water with 0.1% formic acid
0-2	5%	0%
2-5	10%	5%
5-5.5	90%	10%
5.5-8	5%	90%
8-10	5%	90%
10-12	5%	0%

### Mass Spectrometer:

- Positive ion mode electrospray ionization
- Monitor the ion transition of the parent ion ( $m/z$  167) to the product ion ( $m/z$  85) in multiple reaction monitoring (MRM)

Mass Spectrometry Conditions for Cyromazine Quantitation	
capillary voltage	3.1 kV
cone voltage	65 V
collision energy	21-24 V
source temperature	120° C
desolvation temperature	350° C
cone gas flow	135 L/h
desolvation gas flow rate	750 L/h
collision gas	Argon
parent ion	( $m/z$ ) 167
product ion	( $m/z$ ) 85



## References:

- [1] \*Summarized with permission from Xia, Kang, Atkins, Jack et al, "Analysis of Cyromazine in Poultry Feed Using the QuEChERS Method | Coupled with LC-MS/MS" J. Agric. Food Chem, DOI:10.1021/jf9034282
- [2] Listing of instrument manufacturers does not constitute endorsement by UCT

**DCN-012101-196**

UCT, LLC • 2731 Bartram Road • Bristol, PA 19007 800.385.3153 • 215.781.9255

[www.unitedchem.com](http://www.unitedchem.com) Email: [methods@unitedchem.com](mailto:methods@unitedchem.com)

©UCT, LLC 2010 • All rights reserved

