

Extraction of Polycyclic Aromatic Hydrocarbons from Fish Using the QuEChERS Approach



UCT Part Numbers

ECMSSC-MP

4000 mg mag. Sulfate, 1000 mg sodium chloride

ECMPSC1815CT

900 mg MgSO₄, 300 mg PSA and 150 mg endcapped C18

ECPAHR50CT

ENVIRO-CLEAN® PAH centrifuge tubes

Procedure:

1. Extraction: Protect samples from exposure to light to avoid degradation

- 5 grams of homogenized fish are transferred to a 50 mL centrifuge tube.
- Add 10 mL of acetonitrile.
- Mix by shaking.
- Add the contents of pouch containing 4 grams MgSO₄ and 1 gram of NaCl to the centrifuge tube.
- Immediately vortex the mixture for 3 minutes.
- Centrifuge samples for 3 minutes at 3400 rpm.
- Recover the clear supernatant for cleanup.

2. Clean-up, Dispersive Solid-phase (dSPE)

- Use 3 mL of the supernatant for clean-up.
- Add supernatant to centrifuge tube **ECMPSC1815CT**.
- Shake for 1 minute.
- Centrifuge for 1 minute at 3400 rpm.
- Filter supernatant through 0.20 µm PTFE membrane filter
- Samples are ready for analysis.

3. Elution

- Separation of the compounds is performed in a C18 column (CC 150/4 Nucleosil 100-5 C18 PAH, 150 x 4.0 mm; 5 µm particle size; Macherey-Nagel, Duren, Germany) or equivalent maintained at room temperature.
- Inject 15 µL.
- The initial composition of the mobile phase is 50% of ACN and 50% water.
- A linear gradient to 100% is programmed in 15 minutes.
- Final hold 13 minutes.
- Initial conditions are achieved within 1 minute and maintained for 6 min before next run.
- Total run time is 40 minutes.
- Flow rate of 0.8 mL/minute.
- Fluorescence wavelength program:

Each compound is detected at its optimum excitation/emission wavelength pair:

315/260 nm naphthalene, acenaphthene and fluorene

366/260 nm Phenanthrene

430/260 nm anthracene, fluoranthene, pyrene, benz(a)anthracene, chrysene, benzo(b)fluoranthene, benzo(k)fluoranthene, benzo(a)pyrene, dibenz(a,h)anthracene, benzo(g,h,i)perylene and dibenzo(a,l)pyrene

505/290 nm (indeno(1,2,3-cd)pyrene)



References:

- [1] *Adapted from, Ramalhosa, Maria Joao et al, "Analysis of polycyclic aromatic hydrocarbons in fish: evaluation of a quick, easy, cheap, effective, rugged, and safe extraction method", J. Sep. Sci. 2009, 32, 3529 – 3538

DCN-117090-116

UCT, LLC • 2731 Bartram Road • Bristol, PA 19007 800.385.3153 • 215.781.9255

www.unitedchem.com Email: methods@unitedchem.com

©UCT, LLC 2010 • All rights reserved

