


Determination of Carisoprodol and Meprobamate in Whole Blood and Urine Samples by LC/MS-MS



UCT Part Numbers

SLDAGDC21-5UM
SELECTRA® DA (10 mm x 2.1 mm)
Guard Column

SLDA100ID21-5UM
Selectra® DA (100 mm x 2.1 mm,
5 µm) LC Column

CSDAU206
CLEAN SCREEN® DAU
200 mg / 6 mL

SPPHO6001-5
0.1M Phosphate Buffer pH 6

Procedure:

1. Prepare Sample

- To 1 mL of whole blood or urine add internal standard(s)* and 3 mL of 100 mM phosphate buffer (pH 6.0).
- Mix/vortex.
- Sample pH should be 6.0 ± 0.5 .
- Centrifuge for 10 minutes at 3000 rpm.
- Discard pellet.

2. Condition Clean Screen® Extraction Column

- 1 x 3 mL CH₃OH.
- 1 x 3 mL D.I. H₂O.
- 1 x 1 mL 100 mM phosphate buffer (pH 6.0).

Note: Aspirate at < 3 inches Hg to prevent sorbent drying.

3. Apply Sample

- Load sample at 1 to 2 mL/minute

4. Wash Column

- 1 x 3 mL DI H₂O.
- 1 x 1 mL 100 mM acetic acid.
- Dry column (5 minutes at > 10 inches Hg).
- 1 x 2 mL hexane.

5. Elute Carisoprodol/Meprobamate

- 1 x 3 mL hexane/ethyl acetate (50:50);
Collect eluate at 1 to 2 mL / minute.

6. Dry Eluates

- Evaporate to dryness under nitrogen < 40°C
- Reconstitute in 100 µL of mobile phase



Instrument Conditions

Column	100 x 2.1 mm (5 µm) SELECTRA® DA (UCT, Inc.,)
Column Temperature	40 °C
Mobile phase	Isocratic 50:50 CH ₃ OH (0.1% Formic acid): DI H ₂ O (0.1% Formic acid)
Flowrate	0.5 mL/ minute
Detector	API 4000 MS/MS


Compound	Q1	Q3	Time /ms	DP /volts	EP /volts	CXP /volts	CE /volts
Carisoprodol (1)	261.2	55.8	200	16	11	4	47
Carisoprodol (2)	261.2	176.2	200	16	11	4	13
[†] Carisoprodol-d7 (1)	268.2	183.2	200	31	11	4	15
Meprobamate (1)	219.2	158.2	200	26	8.5	4	13
Meprobamate (2)	219.2	54.7	200	26	8.5	4	25
[†] Meprobamate-d7 (1)	226.3	165.2	200	31	9.55	4	13

*Time = Dwell Time; *DP= Declustering Potential; *EP= Exit Potential; *CXP= Collision Cell Exit Potential;

*CE = Collision Energy

[†]Internal Standard

Chromatogram of: Meprobamate, Carisoprodol



DCN-314180-278

UCT, LLC • 2731 Bartram Road • Bristol, PA 19007 800.385.3153 • 215.781.9255

www.unitedchem.com Email: methods@unitedchem.com

©UCT, LLC 2013 • All rights reserved

