


# THC, THC-OH, and Carboxy-THC In Equine Urine for GC/MS Confirmations



## UCT Part Numbers

### XRDAH206

XtrackT DAU  
200 mg, 6 mL

### SPHACE5010-10

Select pH Buffer Pouches  
1M Acetate pH 5.0

### SBSTFA-1-1

SELECTRA-SIL® BSTFA  
w/ 1% TMCS

### BETA-GLUC-10

SELECTRAZYME®  
β – Glucuronidase

## Procedure:

### 1. Prepare Sample

- To 1 mL of urine add internal standard (s) and 50 µL of Beta Glucuronidase solution (*Halotis rufescens*), add 2 mL of 1 M Acetate buffer pH= 5.
- Mix and incubate at 65 °C for 3 hours.
- Cool to room temperature.
- Add 100 of 10 M NaOH. Mix/vortex.
- Hydrolyze for 20 minutes at 60°C.  
Cool before proceeding.
- Adjust sample pH to 3.0 with approx. 1.0 mL of glacial acetic acid. Check pH to insure that the pH value is ~ 3.0. Centrifuge as appropriate.

### 2. Condition XTRACT® DAU Extraction Column

- 1 x 3 mL CH<sub>3</sub>OH.
  - 1 x 3 mL D.I. H<sub>2</sub>O.
  - 1 x 1 mL Acetate buffer (pH= 3.0).
- Note:** Aspirate at full vacuum or pressure.

### 3. APPLY SAMPLE

- Load at 1 to 2 mL/ minute.

### 4. Wash Column

- 1 x 2 mL D.I. H<sub>2</sub>O.
- 1 x 2 mL 100 mM HCl/acetonitrile (95:5).
- Dry column (5-10 minutes at full vacuum or pressure).
- 1 x 200 1 mL hexane; Aspirate.  
(Additional step to remove any residual moisture)

### 5. Elute Cannabinoids

- 1 x 3 mL hexane/ethyl acetate/ glacial acetic acid (49:49:2).
- Collect eluate at 1 to 2 mL/minute.  
**Note:** Before proceeding, ensure there are no water drop lets at the bottom of the collection tube.  
This may increase drying time and decrease BSTFA derivatizing efficiency.


### 6. Dry Eluate

- Evaporate to dryness at < 40 °C.

### 7. Derivatize

- Add 50 µL ethyl acetate and 50 µL BSTFA w/1% TMCS.
- Mix/vortex.
- React 20 minutes at 70°C.
- Remove from heat source to cool.  
**Note:** Do not evaporate BSTFA
- Inject 1 to 2 µL onto gas chromatograph.



**Chromatogram:****THC, THC-OH, and THC-COOH****Mass Spec Table**

| Compound                     | Primary Ion* | Secondary Ion | Tertiary Ion |
|------------------------------|--------------|---------------|--------------|
| THC-TMS                      | 371          | 343           | 366          |
| THC-D3-TMS <sup>†</sup>      | 374          | 346           | 889          |
| THC-OH-TMS                   | 371          | 459           | 474          |
| THC-OH-D3-TMS <sup>†</sup>   | 374          | 462           | 471          |
| THC-COOH-TMS                 | 371          | 473           | 488          |
| THC-COOH-D3-TMS <sup>†</sup> | 374          | 476           | 491          |

\*Quantitation Ion

<sup>†</sup>Suggested internal standard for GC/MS

**4107-03-01**

UCT, LLC • 2731 Bartram Road • Bristol, PA 19007 800.385.3153 • 215.781.9255

[www.unitedchem.com](http://www.unitedchem.com) Email: [methods@unitedchem.com](mailto:methods@unitedchem.com)

©UCT, LLC 2014 • All rights reserved

