


# Extraction and Analysis of Nicotine, Cotinine and Anabasine from Urine using SPE and GC-MS



## UCT Part Numbers

**CSDAU206**  
Clean Screen DAU,  
200mg / 6 mL tube

**SPHPHO6001-5**  
Select pH Buffer Pouch,  
phosphate buffer pH 6

## Procedure:

### 1. Sample Preparation

- Add 1mL of pH 6 phosphate buffer (0.1 M) into a clean sample tube and add internal standards.
- Add 1 mL of urine and mix.
- 2 mL of 0.1 M pH 6 phosphate buffer and mix.
- Centrifuge for 10 minutes at 3000 rpm.

### 2. SPE Column Conditioning

- 1 x 3 mL CH<sub>3</sub>OH.
- 1 x 3 mL D.I. H<sub>2</sub>O.
- 1 x 1 mL 0.1 M phosphate buffer (pH 6.0).

**Note:** Aspirate at < 3 inches Hg to prevent sorbent drying.

### 3. Apply Sample

- Load sample at 1 to 2 mL/minute.

### 4. Wash Column

- 1 x 3 mL D.I. H<sub>2</sub>O.
- 1 x 3 mL 1 M Acetic Acid
- 1 x 3 mL Methanol
- Dry column (5 minutes at > 10 inches Hg or full flow on a positive pressure manifold).

### 5. Elute Nicotine, Cotinine and Anabasine

- Add 1 x 3 mL Methylene chloride: Isopropanol: ammonium hydroxide (78:20:2).

**Note:** Prepare a fresh solution of the Methylene chloride: Isopropanol: ammonium hydroxide mixture by adding 2 mL of ammonium hydroxide to 20 mL of isopropanol. Mix and add 78 mL of methylene chloride.

- Collect eluate at a rate of 1 to 2 mL / minute.

### 6. Dry Eluate

- Evaporate to dryness under nitrogen < 35°C.
- Reconstitute with 100 µL of methanol and vortex mix.
- Transfer to an autosampler vial containing a low volume insert (150 µL)



Instrument Conditions	
Instrumentation	ThermoScientific Trace 1300 GC TriPlus RSH autosampler with a 10 µL syringe
Column	TG-1MS 30 m x 0.25 mm (0.25 µm df)
Carrier Gas	He, Ultrahigh purity
Flow Rate	1.2 mL/ minute
Split Flow	12 mL/ minute
Split Flow Ratio	10:1
Injector Temp	250°C
Septum Purge Flow	5.0 mL/ minute
Transfer Line Temp	280 °C
Injection Volume	1 µL
Oven Run Program	
Initial Temp	50°C
Initial Hold Time	0.5 minute
Temperature Program	30°C / minute
Final Temperature	320°C
Final Hold Time	0.5 minute


Mass Spectrometry	
Detector	Thermo Scientific ISQ
Transfer Line	280°C
Ion Source Temperature	300°C
Ionization Mode	El <sup>+</sup>
Mode	Selective Ion Mode (SIM)
Data Analysis	Xcaliber version 2.2 SPI 2
Note: Quantification was performed on last saved autotune	

Mass Spec Table				
Compound	Retention Time (minute)	Primary Ion	Secondary Ion	Tertiary Ion
Nicotine	5.69	84*	136	162
Nicotine-d4	5.60	136	166*	
Cotinine	7.10	98*	119	176
Cotinine-d3	7.10	101*	122	
Anabasine	6.34	84*	105	133


\*Selected monitoring ion




## Chromatogram of Nicotine (500 ng / mL)



## Chromatogram of Cotinine (500 ng / mL)



## Chromatogram of Anabesine (500 ng / mL)



**4107-01-01**

UCT, LLC • 2731 Bartram Road • Bristol, PA 19007 800.385.3153 • 215.781.9255

[www.unitedchem.com](http://www.unitedchem.com) Email: [methods@unitedchem.com](mailto:methods@unitedchem.com)

©UCT, LLC 2014 • All rights reserved

