

Extraction of Z-Drug/Metabolite Panel From Urine Using Clean Screen® XCEL I SPE and LC-MS/MS



UCT Part Numbers

SLDA100ID21-3UM
Selectra® DA HPLC Column
100 x 2.1, 3µm

SLDAGDC21-3UM
Selectra® DA Guard Column
10 x 2.0, 3µm

Clean Screen® CSXCE106
130 mg/6mL

SLGRDHLDL-HPOPT
Guard Column Holder

Summary:


Z-drugs or nonbenzodiazepines are a class of psychoactive drugs whose pharmacological actions are similar to benzodiazepines and have shown efficacy in treating sleep disorders. Although chemically different, they have similar benefits, risks, and side effects as benzodiazepines. Z-drugs bind to GABA receptor complexes close to or coupled with benzodiazepine receptors. These compounds have been noted as more selective. These compounds bind mainly to the $\alpha 1$ GABA receptor subtype, which mediates hypnotic effects. Benzodiazepines worsen sleep architecture, whereas Z-drugs have little to no impact on sleep architecture, making them preferable.

The three primary groups of Z-drugs are:


- Zolpiclone (Lunesta, Imovane, Zimovane, Imrest)
- Zolpidem (Ambien, Ambien CR, Stilnoct, Intermezzo, Stilnox)
- Zaleplon (Sonata, Starnoc)

Zolpidem is partially metabolized in the liver to form an inactive N-demethylated derivative and an active N-oxide metabolite. These metabolites account for about 30% of the initial dose in urine, while 7-10% of the drug is excreted unchanged. Less than 1% of Zolpidem and Zaleplon is excreted in urine unchanged. Due to these low concentrations, sample cleanup using SPE is usually required to eliminate undesirable matrix and concentrate the sample.


Zolpidem
pKa=5.65



Zolpiclone
pKa=6.89



Zaleplon
pKa=0.3



Clean Screen® XCEL I columns were chosen because of their excellent capability to extract basic compounds while eliminating the need for time-consuming column conditioning. This method reduces the time needed to extract a panel of samples and minimizes solvent use for sample cleanup. To ensure the functional groups on the various compounds were fully ionized, the urine samples were adjusted to a pH of 4 before passing them through the column. At this pH, samples were effectively retained on the column. A wash step was performed, which removed unwanted matrix and interferences without losing any of the analytes in question. The Z-drugs were eluted using a highly basic solvent combination of MeCl₂: IPA: NH₄OH (pH 11-12).

Most of the compounds in the evaluated panel are true bases with ionizable functional groups. However, Zaleplon is structurally more similar to a benzodiazepine. Because this drug only becomes charged in extremely acidic pH conditions, it functions more as a neutral compound. Retained primarily via hydrophobic interactions, the amount of organic in the wash needed to be optimized not to compromise the overall recovery of Zaleplon but still provide sufficient cleanup for the remainder of the analytes.

Procedure

1. Sample Preparation

- To 1 mL of urine sample, add 3 mL of 0.1% Formic Acid Solution.
- Vortex Sample.

2. Apply Sample to Clean Screen® XCEL I column

- Load sample directly to column without any preconditioning.
- Pull sample through at a rate of 1-2 mL/ minute.
- Dry column thoroughly under full vacuum or positive pressure for 1 minute.

3. Wash

- 1 x 3 mL 75:25 100mM Acetic Acid:MeOH.
- 1 x 3 mL Hexane.
- Dry column thoroughly under full vacuum or positive pressure for 5-10 minutes.

4. Elution

- 1 x 3 mL CH₂Cl₂/IPA/NH₄OH (78:20:2).
- Collect eluate at 1 to 2 mL/minute.

Note: Prepare elution solvent daily. Add IPA/NH₄OH, mix, then add CH₂Cl₂ (pH 11-12).

5. Dry Eluate

- Evaporate the fraction to complete dryness under stream of dry air or nitrogen at ~ 35 °C.

6. Reconstitute

- Reconstitute the sample in 100 µL of mobile phase.



Instrumental:


HPLC Conditions	
Instrumentation	Agilent 1200 Series Binary Pump SL
HPLC column	Selectra® DA, 100 × 2.1 mm, 3 µm
Guard column	Selectra® DA, 10 × 2.0 mm, 3 µm
Guard column holder	p/n: SLDGRDHLDR
Column temp.	50°C
Autosampler temp.	10°C
Injection volume	10 µL
Mobile phase A	0.1% formic acid in H ₂ O
Mobile phase B	0.1% formic acid in ACN
Flow rate	300 µL/min

LC Gradient		
Time (min)	A (%)	B (%)
0	95	5
0.5	80	20
6	80	20
7.5	0	100
12.0	0	100
13.2	95	5
18	95	5

MS Conditions	
Instrumentation	API 4000 QTRAP MS/MS
Ionization mode	ESI ⁺
Spray voltage	4200 V
Vaporizer temperature	650 °C
Sheath gas pressure	40 arbitrary units
Auxiliary gas pressure	5 arbitrary units
Collision gas pressure	Medium



Chromatogram:



Analyte		MRM Transitions		RRT (min)
		Q1	Q3	
1	N-Desmethylopiclone	375.033	244.900	7.36
2	Zopiclone	388.988	244.900	7.86
3	Zopiclone-N-oxide	404.966	143.000	9.50
4	Zolpidem	308.289	234.900	10.12
5	Zaleplon	305.963	264.000	10.80

Results:

Compound	50 ng/mL (n=5)		300 ng/mL (n=5)	
	Extraction Recovery	Matrix Effect	Extraction Recovery	Matrix Effect
N-Desmethyl Zopiclone	88%	8%	89%	5%
Zopiclone	90%	44%	91%	23%
Zopiclone-N-Oxide	71%	46%	75%	37%
Zolpidem	80%	26%	93%	26%
Zaleplon	90%	13%	102%	20%



Tip:

[1] *Insignificant loss of Zaleplon was observed using 3mL 75:25 100mM Acetic Acid:MeOH wash in this study. To correct for any residual recovery or matrix issues, the use of deuterated internal standards is strongly recommended.

5104-03-02

UCT, LLC • 2731 Bartram Road • Bristol, PA 19007 800.385.3153 • 215.781.9255

www.unitedchem.com Email: methods@unitedchem.com

©UCT, LLC 2015 • All rights reserved

