


# Determination of Polycyclic Aromatic Hydrocarbons in Drinking Water by Liquid-Solid Extraction and HPLC with Coupled Ultraviolet and Fluorescence Detection



## UCT Part Numbers

<b>ECUNIPAH</b> 2000 mg unendcapped C18, 83 mL cartridge	Or	<b>EUC1812M15</b> 2000 mg unendcapped C18, 15 mL cartridge
<b>CUC181M6</b> 1000 mg unendcapped C18, 6 mL cartridge		<b>ECSS25K</b> Anhydrous Sodium Sulfate

## Procedure:

### 1. Cartridge Preparation

- Wash with 4 x 10 mL aliquots of methylene chloride ( $\text{MeCl}_2$ ).
- Wash with 4 x 10 mL aliquots of methanol ( $\text{MeOH}$ ).
- Wash with 2 x 10 mL aliquots of reagent water.

**Note:** Do not let the cartridge dry out after step 1) c. otherwise repeat starting at 1) b.

### 2. Sample Extraction

- Adjust the vacuum setting for a flow rate of 10-15 mL per minute.
- Add the 1 liter sample to the cartridge.
- Rinse sample bottle with reagent water, add to cartridge and draw through.
- Dry cartridge by drawing full vacuum for 10 minutes.

### 3. Sample Elution and Drying


- Elute the cartridge dropwise by using 2 x 5 mL aliquots of  $\text{MeCl}_2$  and collect.
- Rinse sample container with 5 mL of  $\text{MeCl}_2$ , add to cartridge and draw through.
- Prepare a drying column/funnel containing 10-20 g sodium sulfate by rinsing with 10 mL of  $\text{MeCl}_2$  and discard.
- Add the eluate to the drying column, draw through and collect.
- Rinse the eluate vial and drying column with a 2 x 5 mL aliquot of  $\text{MeCl}_2$  and collect.

### 4. Sample Evaporation

- Evaporate the extract using a gentle stream of  $\text{N}_2$  with a water bath or heating block temperature of 40°C. Evaporate to about 1.0 mL.
- Add 3.0 mL of acetonitrile (ACN).
- Concentrate to a final volume of 0.5 mL.

### 5. Sample Analysis

- Inject 5 - 100  $\mu\text{L}$  into an HPLC.



## References:

[1] See "Determination of Polycyclic Aromatic Hydrocarbons in Drinking Water by Liquid-Solid Extraction and HPLC with Coupled Ultraviolet and Fluorescence Detection", W. J. Bashe & T.V. Baker (Technology Applications, Inc, Environmental Monitoring Systems Laboratory, US Environmental Protection Agency, Cincinnati, OH)

**DCN-017009-191**

UCT, LLC • 2731 Bartram Road • Bristol, PA 19007 • 800.385.3153 • 215.781.9255

[www.unitedchem.com](http://www.unitedchem.com) Email: [methods@unitedchem.com](mailto:methods@unitedchem.com)

©UCT, LLC 2012 • All rights reserved

