


Determination of Nitrosamines in Drinking Water by Solid Phase Extraction and Capillary Column Gas Chromatography with Large Volume Injection and Chemical Ionization Tandem Mass Spectrometry (MS/MS)*



UCT Part Numbers

EU52112M6

2000 mg activated coconut carbon, 6 mL

ECSS156


6 mL Drying Cartridge with 5 grams anhydrous sodium sulfate

Summary:

Activated carbon is used for the determination of various nitrosamines in finished drinking water and untreated source waters using GC/MS/MS.

Analyte	Abbreviation	% Recovery n=3
N-Nitrosodimethylamine	NDMA	95
N-Nitrosomethyldiethylamine	NMEA	98
N-Nitrosodiethylamine	NDEA	95
N-Nitrosodi-n-propylamine	NDPA	90
N-Nitrosodi-n-butylamine	NDBA	94
N-Nitrosopyrrolidine	NPYR	76
N-Nitrosopiperidine	NPIP	81

*For complete details on Method 521, September 2004, the analyst is referred to: J.W.Munch & M.V.Bassett, "Determination of Nitrosamines in Drinking Water by Solid Phase Extraction and Capillary Column Gas Chromatography with Large Volume Injection and Chemical Ionization Tandem Mass Spectrometry (MS/MS), National Exposure Research Laboratory, Office of Research and Development, U.S. Environmental Protection Agency, Cincinnati, OH 45268



Procedure:

1. Cartridge Conditioning

- Add 3 mL of methylene chloride to the cartridge, then slowly draw all solvent through the cartridge.
- Add 3 mL of methanol to the cartridge, turn on vacuum and draw through.
- Add 3 mL of methanol again and draw through so that the methanol just covers the top of the cartridge frit.

Note: Do not let the cartridge go dry after this step otherwise repeat starting at step 1 b)

- Add 3 mL of reagent water and draw through.
- Repeat water rinse, step d) **5 additional times**.

Note: Proper conditioning of the cartridge is essential for good precision and accuracy

2. Sample Extraction

- Adjust the vacuum setting so that the flow rate is 10 mL/minute.
- After sample extraction, draw air through the cartridge for 10 minutes at full vacuum.
- After drying, proceed immediately to cartridge elution


3. Cartridge Elution

- Insert a clean collection tube in the manifold.
- Fill the cartridge with methylene chloride.
- Partially draw the methylene chloride through at low vacuum and then turn vacuum off.
- Allow cartridge to soak for 1 minute.
- Draw the remaining methylene chloride through in dropwise fashion.
- Continue to add methylene chloride to the cartridge as it is being drawn through until a total of 12-13 mL have been added.
- Concentrate the methylene chloride to about 0.9 mL in a water bath near room temperature. Do not concentrate less than 0.5 mL as loss of analyte may occur.

Note: Small amounts of residual water from the sample container and SPE cartridge may form an immiscible layer with the extract. To eliminate the water a drying column packed with 5 grams of anhydrous sodium sulfate or use ECSS15M6 for drying. Wet the cartridge with a small volume of methylene chloride before adding extract. Rinse the drying column with 3 mL of methylene chloride.

4. Sample Analysis

- Calibrate the MS in EI mode using FC-43.
- Inject into a GC/MS/MS.
- Identify the product ion spectrum to a reference spectrum in a user created database.



Mass Spectral Data

Analyte	Retention Time (min)	Precursor Ion (m/z)	Product/Quantitation Ion (m/z)
NDMA	8.43	75	43(56)
NMEA	11.76	89	61(61)
NDEA	14.80	103	75(75)
NPYR	22.34	101	55(55)
NDPA	22.40	131	89(89)
NPIP	24.25	115	69(69)
NDBA	30.09	159	57(103)
NDMA-d6 surrogate	8.34	81	46(59)
NMEA-d10 IS	14.63	113	81(81)
NDPA-d6 IS	22.07	145	97(97)

Injector Program


Temp (°C)	Rate (°C/min)	Time (min)
37	0	0.72
250	100	2.13
250	0	40

Injector Split Vent Program

Time (min)	Split Status
0	Open
0.70	Closed
2.00	Open

GC Oven Temperature Program

Temperature (°C)	Rate (°C/min)	Hold Time (min)
40	0	3.0
170	4.0	0



Limits and Lowest Concentration Minimum Reporting Levels

Analyte	DL (ng/L)	LCMRL (ng/L)
NDMA	0.28	1.6
NMEA	0.28	1.5
NDEA	0.26	2.1
NPYR	0.35	1.4
NDPA	0.32	1.2
NPIP	0.66	1.4
NDBA	0.36	1.4

DCN-216111-112

UCT, LLC • 2731 Bartram Road • Bristol, PA 19007 800.385.3153 • 215.781.9255

www.unitedchem.com Email: methods@unitedchem.com

©UCT, LLC 2012 • All rights reserved

