

Determination of Ephedra Alkaloids & Synephrine in Dietary Supplements via Strong Cation-Exchange SPE and LC-MS/MS Detection



UCT Part Numbers

CUBCX1HL56

High-load Benzenesulfonic Acid
500 mg / 6 mL SPE column

SLPFPPGDC20-3UM

Selectra® PFPP guard cartridge
10 × 2.1 mm, 3 μm

SLPFPP100ID21-3UM

Selectra® PFPP HPLC column
100 × 2.1 mm, 3 μm

SLGRDHLDR-HPOPT

Guard cartridge holder

Summary:

Ephedra alkaloids are phenethylamines that occur naturally in plants, including the herb Ma Huang used in traditional Chinese medicine. Ephedra alkaloids are potent CNS stimulants and also have a sympathomimetic effect on the peripheral nervous system. Some active ingredients of these plants are used as ingredients in cold remedies (e.g. pseudoephedrine). The ephedra alkaloids are also incorporated into dietary supplements to promote weight loss or to increase alertness and physical activity (e.g. body building). However, severe contraindications have been reported for individuals with hypertension or other cardiovascular diseases, particularly when used in combination with caffeine [1]. Products containing ephedrine were popular dietary supplements until the FDA banned their use in 2004 [2]. Since then the active ingredient in dietary supplements has largely been replaced by synephrine, a naturally occurring alkaloid found in plants such as Citrus fruits; it is similar in structure to ephedrine.

The ephedra alkaloids are small, hydrophilic, basic analytes that are difficult to retain and separate on traditional HPLC columns using alkyl-bonded stationary phases. They are capable of strongly interacting with free silanols on the surface, which leads to peak tailing of the analytes and affects the resolution and quantification.

Current methodologies used for the separation of ephedra alkaloids use reversed-phase columns with ion-pairing reagents, time-consuming derivatization procedures, or use strong cation-exchange phases. However, these approaches are not very amenable to LC-MS/MS analysis. An alternative approach to traditional alkyl phases is the use of a fluorinated stationary phase. In addition to dispersive interactions available on traditional alkyl phases, pentafluorophenylpropyl phases can undergo dipole-dipole, and pi-pi interactions. This imparts unique selectivity to the column that can sufficiently resolve the ephedra alkaloids.

The aim of this study was to develop a multi-analyte procedure for the extraction, cleanup, and quantification of the ephedra alkaloids in functional foods and natural products. High capacity strong cation-exchange SPE cartridges were used for the isolation of the phenethylamines from dietary supplements. HPLC separation, including separation of the stereoisomers, was carried out using a UCT Selectra® PFPP column prior to detection by LC-MS/MS.



Sample Pretreatment:

Weigh 1 ± 0.1 g of dietary supplement in question into a 15 mL polypropylene centrifuge tube. For this study, a SPEX® SamplePrep® GenoGrinder® was used to pulverize the tablets. Add 10 mL of 1% formic acid to each sample. Shake or vortex sample for 15 minutes to fully extract the ephedra alkaloids. Ensure samples are fully dissolved. Centrifuge the sample for 10 min at $\geq 3000 \times g$ and 4°C .

SPE Procedure:

1. SPE Conditioning

- Add 2×4 mL of methanol to CUBCX1HL56 SPE cartridge.
- Add 4 mL of ultrapure water.
- Add 4 mL of 1% formic acid

Note: Do not let the cartridge go dry otherwise repeat steps a) through c).

2. Sample Extraction

- Load supernatant from step 1d).
- Allow sample to percolate through the cartridge or apply a vacuum if necessary (adjust vacuum for flow of 1–3 mL per minute).

3. Second Extraction (Optional)

- Add 5 mL of 1% formic acid to each sample.
- Shake or vortex sample for 5 minutes.
- Centrifuge the sample for 10 min at $\geq 3000 \times g$ and 4°C .
- Apply supernatant to the SPE cartridge.

4. Wash Cartridge

- Add 2×4 mL of 0.1% formic acid and slowly draw through.
- Add 2×4 mL methanol and slowly draw through.
- Dry under vacuum for ≈ 30 sec to remove excess solvent.

5. Elution

- Elute the ephedra alkaloids using 8 mL of methanol containing 2% ammonium hydroxide.
- Evaporate off the methanol solvent at 40°C under a gentle stream of nitrogen until it reaches a volume of ≈ 1 mL.
- Add 1 mL of aqueous mobile phase (10mM ammonium acetate).
- Evaporate off any remaining methanol.

Note: Ephedra alkaloids are similar to amphetamines, which are known to be volatile compounds. Therefore, extra care was taken during the evaporation step to avoid any potential loss in recovery that may occur during this step.

- Vortex the samples for 1 min and filter through a $0.2 \mu\text{M}$ syringe filter directly into a HPLC vial.



LC-MS/MS Parameters:

Table 1. Instrumentation	
HPLC System	Thermo Scientific™ Dionex™ Ultimate™ 3000
MS System	Thermo Scientific™ TSQ Vantage™ (MS/MS; APCI ⁺)
Ionization Mode	ESI ⁺
HPLC Column	UCT Selectra® PFPP, 100 × 2.1 mm, 3 μm (p/n: SLPFPP100ID21-3UM)
Guard Column	UCT Selectra® PFPP, 10 × 2.1 mm, 3 μm, (p/n: SLPFPPGDC20-3UM)
Guard Column Holder	p/n: SLGRDHLDR
Column Temperature	50 °C
Flow Rate	500 μL/min
Injection Volume	5 μL
Isocratic Elution	85:15 (A:B, v:v)

Table 2. MRM Transitions							
Compound	t _R (min)	Precursor ion	Product ion 1	CE 1	Product ion 2	CE 2	S-lens (V)
Ephedrine	9.9	166.17	91.01	29	115.02	22	46
Pseudoephedrine	11.4	166.17	91.01	29	115.02	22	46
Norephedrine	5.6	152.141	91.02	31	115.02	22	38
Norpseudoephedrine	6.7	152.141	91.02	31	115.02	22	38
Methylephedrine	14.7	180.17	90.99	31	147.04	15	52
Synephrine	2.1	168.10	90.97	19	106.97	29	37
Ephedrine-d ₃ (IS)	9.9	169.17	90.99	31	115.00	24	49
Pseudoephedrine-d ₃ (IS)	11.4	169.17	90.99	31	115.00	24	49

Results and Discussion:

Table 3. Accuracy & precision Data at 100 ppb (n=5)						
	Ephedrine	Pseudoephedrine	Norephedrine	Norpseudoephedrine	Methylephedrine	Synephrine*
Sample 1	94.13	93.07	80.22	108.30	79.83	47.35
Sample 2	91.62	94.28	48.49	88.11	79.91	38.14
Sample 3	92.62	92.74	63.24	89.02	89.72	52.41
Sample 4	92.60	93.56	61.63	75.05	79.48	45.43
Sample 5	93.39	93.56	49.86	70.75	83.24	51.85
Mean	92.87	93.44	60.69	86.25	82.44	47.04
RSD	1.02	0.62	21.08	17.03	5.27	12.30

*Note: Synephrine is more polar than the ephedra alkaloids and is not as well retained on the sorbent. It is recommended to include an isotopically labeled internal standard in order to achieve the best results.






Figure 1. Example Calibration Curve (ephedrine) over a 10-250 ng/mL concentration range.




Figure 2. Chromatogram of a 10 ng/mL standard



Conclusions:

A method was successfully developed for the extraction, cleanup, and quantification of the ephedra alkaloids and synephrine in functional foods and natural products. Strong cation-exchange SPE was used to isolate the phenethylamines from the complex sample matrix, which consisted of 9 herbal extracts containing a high concentration of calcium, caffeine and additional excipients. A high capacity SPE sorbent was used as it offers better retention than standard SCX sorbent. It is recommended to include isotopically labeled internal standards into the method, particularly for the hydrophilic synephrine which is not as efficiently retained by SPE as the ephedra alkaloids. HPLC separation of the 5 ephedra alkaloids and synephrine was successfully conducted on UCT's Selectra® PFPP column, including baseline resolution of the 2 sets of stereoisomers included in the method (ephedrine / pseudoephedrine and norephedrine / norpseudoephedrine). Good LC-MS/MS sensitivity was observed for all compounds (< 10 ng/mL).

References:

[1] Journal of Chromatography A, 1161 (2007) 71–88

[2] Food and Drug Administration, Federal Registry 69 (2004) 6787

7106-01-03

UCT, LLC • 2731 Bartram Road • Bristol, PA 19007 800.385.3153 • 215.781.9255

www.unitedchem.com Email: methods@unitedchem.com

©UCT, LLC 2017 • All rights reserved

