



Extraction and Analysis of Nicotine, Cotinine and Anabasine from Urine using SPE and GC-MS

UCT Part Numbers:

CSDAU206: Clean Screen DAU, 200mg / 6 mL tube

SPPHO6001-5: Select pH Buffer Pouch, phosphate buffer pH 6

July 2014

Procedure

1. Sample Preparation

- a) Add 1mL of pH 6 phosphate buffer (0.1 M) into a clean sample tube and add internal standards
- b) Add 1 mL of urine and mix
- c) 2 mL of 0.1 M pH 6 phosphate buffer and mix
- d) Centrifuge for 10 minutes at 3000 rpm

2. SPE Column Conditioning

- a) 1 x 3 mL CH₃OH.
- b) 1 x 3 mL D.I. H₂O
- c) 1 x 1 mL 0.1 M phosphate buffer (pH 6.0).

NOTE: Aspirate at < 3 inches Hg to prevent sorbent drying.

3. Apply Sample

- a) Load sample at 1 to 2 mL/minute.

4. Wash Column

- a) 1 x 3 mL D.I. H₂O
- b) 1 x 3 mL 1 M Acetic Acid
- c) 1x 3 mL Methanol
- d) Dry column (5 minutes at > 10 inches Hg or full flow on a positive pressure manifold).

5. Elute Nicotine, Cotinine and Anabasine

- a) Add 1 x 3 mL Methylene chloride: Isopropanol: ammonium hydroxide (78:20:2).

NOTE: Prepare a fresh solution of the Methylene chloride: Isopropanol: ammonium hydroxide mixture by adding 2 mL of ammonium hydroxide to 20 mL of isopropanol. Mix and add 78 mL of methylene chloride.

- b) Collect eluate at a rate of 1 to 2 mL / minute.

6. Dry Eluate

- a) Evaporate to dryness under nitrogen < 35°C.
- b) Reconstitute with 100 µL of methanol and vortex mix.
- c) Transfer to an autosampler vial containing a low volume insert (150 µL)

Instrument Conditions

| | | |
|--------------------|--|-------------|
| Instrumentation | ThermoScientific Trace 1300 GC TriPlus RSH autosampler with a 10 µL syringe | |
| Column | TG-1MS 30m x 0.25mm (0.25µm df) | |
| Carrier Gas | He, Ultrahigh purity | |
| Flow Rate | 1.2mL/ minute | |
| Split Flow | 12 mL/ minute | |
| Split Flow Ratio | 10:1 | |
| Injector Temp | 250°C | |
| Septum Purge flow | 5.0 mL/ minute | |
| Transfer Line Temp | 280°C | |
| Injection Volume | 1uL | |
| Oven Run Program | | |
| | Initial Temp | 50°C |
| | Initial Hold Time | 0.5minute |
| | Temperature Program | 30°C/minute |
| | Final Temperature | 320°C |
| | Final Hold Time | 0.5minute |

Mass Spectrometry

| | |
|---|----------------------------|
| Detector | Thermo Scientific ISQ |
| Transfer Line | 280°C |
| Ion Source Temperature | 300°C |
| Ionization Mode | EI+ |
| Mode | Selective Ion Mode (SIM) |
| Data Analysis | Xcaliber version 2.2 SPI 2 |
| Note: Quantification was performed on last saved autotune | |


Mass Spec Table

| Compound | Retention Time (minute) | Primary Ion | Secondary Ion | Tertiary Ion |
|--------------|-------------------------|-------------|---------------|--------------|
| Nicotine: | 5.69 | 84* | 136 | 162 |
| Nicotine-d4: | 5.60 | 136 | 166* | |
| Cotinine: | 7.10 | 98* | 119 | 176 |
| Cotinine-d3: | 7.10 | 101* | 122 | |
| Anabesine: | 6.34 | 84* | 105 | 133 |

*Selected monitoring ion

Chromatogram of Nicotine (500 ng / mL)

RT: 5.00 - 6.69 SM: 15G




NL: 3.43E7

Base Peak m/z= 83.50-84.50 F:
{0,0} + cEI SIM ms
[83.50-84.50, 135.50-136.50,
161.50-162.50] MS ICIS
05072014_Nicotine_Sample3

Chromatogram of Cotinine (500 ng / mL)


RT: 6.73 - 7.48 SM: 15G



NL: 2.18E7
Base Peak m/z= 97.50-98.50 F:
{0,0} + c EI SIM ms
[97.50-98.50, 118.50-119.50,
175.50-176.50] MS
05062014_Nicotine_500ng_m
L

Chromatogram of Anabasin (500 ng / mL)

RT: 5.97 - 6.72 SM: 15G



NL: 4.19E6
Base Peak m/z= 83.50-84.50 F:
{0,0} + c EI SIM ms
[83.50-84.50, 104.50-105.50,
132.50-133.50] MS ICIS
05052014_Nicotine_500ng_mL

4107-01-01