

A Simple SPE Method for the Determination of Malachite Green, Crystal Violet and Other Synthetic Dyes in Seafood

UCT Part Numbers

CSDAU206

Clean Screen[®] DAU 200 mg/6mL Column

SLC-18100ID21-3UM

Selectra[®] C18 HPLC Column 100 x 2.1 mm, 3 μm

SLC-18GDC20-3UM

Selectra[®] C18 HPLC Guard Column 10 x 2.0 mm, 3 μm

SLGRDHLDR

Selectra[®] HPLC Guard Column Holder





Summary:

Triphenylmethane dyes (TPMs) are synthetic dyes used for a wide range of industrial applications. However, they are also used in aquaculture production due to their antibacterial, anti-fungal and anti-parasitic properties. They are cheap, very effective and readily available, but due to their toxicity to humans (mutagenic & carcinogenic) the TPMs are banned in the US and EU. In both regions there is a zero-tolerance policy regarding their use. As such, the respective authorities (FDA and EU Commission) have established very low detection limits for these compounds in seafood products, namely 1 μ g/kg (US) and 2 μ g/kg (EU). As the parent dyes are rapidly metabolized to form their leuco- metabolites that predominate and persist in fish tissue, the detection limits are expressed as the sum of parent drug and corresponding leuco-form.

Analyzing TPMs can be difficult due to their instability, which causes problems during sample preparation and instrumental analysis. The TPM dyes readily undergo oxidation-reduction reactions and are prone to UV degradation (demethylation). As such, it is important to take these issues into account when developing a method or when doing routine analysis of these compounds. To stabilize the dyes, it is common to add an antioxidant prior to sample extraction and in the final sample extract prior to detection. Common antioxidants used include TMPD, hydroxylamine and ascorbic acid. Analytical standards solutions also experience stability issues and should be prepared in an aprotic solvent (e.g. acetonitrile) as the leuco metabolites are not stable for a long period of time in H2O and MeOH (even at -20 °C) and oxidize to the non-leuco form. This can in turn affect the reliability of the analytical method by leading to poor recoveries/reproducibility and non-linear calibration curves.

Most reported methods for TPM analysis use lengthy sample preparation procedures or do little to no sample cleanup. Most methods include malachite green, leucomalachite green, crystal violet, leucocrystal violet and sometimes the structurally related dye brilliant green. HPLC-UV is often used for their analysis due to the compounds good UV absorbance, but this technique can only detect the parent compounds and not the leucometabolites. To include these, an oxidation step needs to be incorporated into the LC system using a post-column reactor to convert the leuco- form to the chromophorous parent form. Alternatively, the use of LC-MS/MS avoids this requirement and is capable of detecting all TPM dyes. It also offers the advantage of simpler sample preparation methodology and the inclusion of additional compounds besides the classical TPM dyes.

Sample preparation approaches for TPM analysis often use time consuming, complicated procedures consisting of numerous steps and requiring large volumes of solvent to evaporate. SPE using strong-cation exchange sorbent is also used and is very effective at removing matrix interferences. However, it is difficult to elute dyes from SCX sorbent as most dyes carry a permanent positive charge (see Figure 1). Most methods using SCX SPE use a high salt concentration (> 1 M) for elution. However, this is not very suitable for LC-MS/MS analysis due to the lack of volatility of such high salt concentrations.

This application note uses an alternative elution solvent composed of methanol (MeOH) containing 1% triethylamine (TEA) and 0.5% formic acid (i.e. triethylammonium formate). The positively charged TEA (a quaternary amine) then acts as a counter-ion to release the dyes from the sorbent. No evaporation step is included in the method which avoids any loss that may occur with this step. The eluted extracts can be analyzed directly as TEA is LC-MS/MS compatible (due to its volatility) and has no ion-pairing effect on the dyes as they carry the same charge. Ascorbic acid is used as antioxidant because positively charged hydroxylamine (quaternary amine) will interfere with the retention on the sorbent (same mechanism as elution solvent). The simple method outlined here allows for the rapid analysis of dyes in seafood while achieving good accuracy/precision and low sensitivity.




Figure 1. Structures of nine compounds included in method

5

Procedure:

1. Sample Preparation

- a) Weigh 2 ± 0.1 g of sample into a 15mL polypropylene centrifuge tube.
- b) Add 10 mL of 1% formic acid in acetonitrile and 1mL of 1M ascorbic acid (antioxidant) to each sample.
- c) Shake on a SPEX[®] SamplePrep[®] GenoGrinder[®] (or alternative mechanical shaker) or vortex samples for 15 minutes to fully extract the dyes.
- d) Centrifuge the samples for 10 min at \geq 3000 rcf and 4°C.
- e) Transfer the supernatant to a 50mL polypropylene centrifuge tube and add 20mL of Mcilvaine's buffer (0.1M, pH 3.5)**.
- f) Vortex the samples briefly (1 min) to mix the organic solvent with the buffer.
- g) Prior to the SPE step, centrifuge the samples for 5 min at ≥3000 rcf and 4°C. This makes the SPE step easier by separating out any solid particles and/or lipids in the sample.

2. Condition Cartridge

- a) Add 3 mL of methanol to SPE cartridge CSDAU256.
- b) Add 3 mL of ultrapure water.
- c) Add 1 mL of Mcilvaine's buffer (0.1M, pH 3.5).

Note: Do not let the cartridge go dry otherwise repeat steps a) through c).

3. SPE Extraction

- a) Load supernatant from step 1
- b) Adjust vacuum for flow of 1–3 mL per minute.

4. Wash Cartridge

- a) Add 3 mL of 0.1% Formic Acid and slowly draw through.
- b) Add 3 mL of 0.1% Formic Acid in MeOH and slowly draw through.
- c) Dry under vacuum for 1 minute to remove excess solvent.

Note: A portion of the dyes will stick to the SPE frit due to the hydrocarbon backbone. The wash steps included in the method uses acidified solutions to maintain a positive charge on the dyes, which ensures they are retained on the cation-exchange sorbent.

5. Elute Cartridge

FOOD

- a) Elute the dyes from the SPE cartridge using 3-4 mL elution solvent (1% TEA + 0.5% formic acid in MeOH).
- b) Vortex the samples for 2 min and transfer a 1 mL aliquot to an autosampler vial for analysis.

** Mcilvaines buffer pH 3.5 – mix 500mL of 0.1M citric acid with 400mL 0.1M disodium hydrogen phosphate. Adjust pH to the correct value (3.5) using the citrate or phosphate solution.

LC-MS/MS Conditions:

	HPLC Conditions							
Instrumentation	Thermo ScientificTM DionexTM UltimateTM 3000 LC system							
HPLC column	electra [®] C18 HPLC Column 100 x 2.1 mm, 3 μm							
Guard column	electra [®] C18 HPLC Guard Column 10 x 2.0 mm, 3 μm							
Run time	min (including re-equilibration time)							
Column temp.	40°C							
Injection volume	5 μL							
Autosampler	10°C							
Wash solvent	MeOH:ultrapure water (1:1, v/v)							
Flow rate	400 μL/min							
Divert valve	Mobile phase was sent to waste (0–3.5 & and 8–12min) to prevent ion source							
Mobile Phase A	0.1% Formic Acid in ultrapure water							
Mobile Phase B	0.1% Formic Acid in MeOH							

Mobile Phase Gradient						
Time	A (%)	B (%)				
0.0	90	10				
2.0	90	10				
3.0	40	60				
4.0	40	60				
5.0	0	100				
8.0	0	100				
8.2	90	10				
12.0	90	10				

	MS Conditions
Instrumentation	Thermo Scientific [™] TSQ Vantage [™] tandem mass spectrometer
Ionization mode	ESI+
Spray voltage	4000 V
Vaporizer temperature	350°C
Capillary temperature	300°C
Sheath gas pressure	60 arbitrary units
Auxiliary gas pressure	5 arbitrary units
lon sweep gas	0 arbitrary units
Declustering potential	0 V
Q1 and Q3 peak width	0.2 and 0.7 Da
Collision gas	Argon
Collision gas pressure	2.0 mTorr
Acquisition method	EZ method (SRM)
Cycle time	0.6 sec
Software for data processing	TraceFinderTM version 3.0
Weighting factor applied to calibration	1/X



SRM Transitions									
Analyte	t _R (min)	Precursor ion	Product ion 1	CE 1	Product ion 2	CE 2	S-lens (V)		
Azure B	4.74	270.1	212.0	48	254.1	31	100		
Methylene Blue	4.77	284.1	252.0	53	268.1	33	100		
Crystal Violet	4.93	374.2	239.0	33	358.2	28	100		
Malachite Green	5.59	329.1	165.1	59	313.2	37	100		
Nile Blue	5.76	318.1	245.1	55	274.1	35	100		
Leucomalachite Green	5.98	331.1	223.1	50	239.1	29	100		
Leucocrystal Violet	5.99	372.2	340.2	51	356.2	37	100		
Brilliant Green	6.19	385.2	297.1	51	341.2	37	100		
Victoria Blue	6.27	470.2	333.1	48	349.1	34	145		

Results:

Accuracy & Precision Data for 9 Dyes Spiked in Salmon Tissue (n=5)									
	MG	LMG	CV	LCV	NB	AB	MB	BG	VB
1 ppb									
Sample 1	92.2	83.8	88.1	101.7	87.2	94.2	76.5	98.5	98.2
Sample 2	94.0	109.5	83.4	92.8	85.4	97.3	80.0	92.9	90.4
Sample 3	92.6	91.2	83.0	95.5	85.0	96.8	77.6	92.1	91.4
Sample 4	93.5	105.6	94.9	91.3	85.1	94.3	73.3	93.2	92.5
Sample 5	92.4	81.2	83.1	103.0	87.9	96.0	76.8	99.0	93.4
Mean	92.9	94.3	86.5	96.9	86.1	95.7	76.9	95.1	93.2
RSD	0.83	13.52	5.97	5.44	1.58	1.51	3.13	3.51	3.26
	MG	LMG	CV	LCV	NB	AB	MB	BG	VB
				10 pp	b				
Sample 1	97.1	91.1	92.6	100.4	95.9	96.6	85.6	106.9	104.4
Sample 2	99.7	96.6	99.1	105.9	98.6	98.6	81.8	103.2	104.4
Sample 3	98.7	94.9	98.3	97.6	97.5	99.4	87.0	104.2	102.7
Sample 4	99.9	88.1	91.0	107.3	99.1	97.9	88.2	98.5	92.2
Sample 5	100.0	97.5	98.6	97.4	99.5	93.5	86.5	102.1	104.3
Mean	99.1	93.7	95.9	101.7	98.1	97.2	85.8	103.0	101.6
RSD	1.24	4.24	3.99	4.56	1.49	2.39	2.85	2.98	5.21





8106-01-01 UCT, LLC • 2731 Bartram Road • Bristol, PA 19007 800.385.3153 • 215.781.9255 www.unitedchem.com Email: methods@unitedchem.com ©UCT, LLC 2018 • All rights reserved