

Analysis of Low-Level Acidic, Basic, & Neutral Drugs of Abuse in Oral Fluids via LC-MS/MS Detection



UCT Part Numbers

CSXCE103

Clean Screen® XCEL I
130 mg, 3 mL Column

SPPHO6001-10

Select PH Buffer Pouches
100 mM Phosphate Buffer pH 6.0

SLFPPP100ID21-3UM

Selectra® PFPP HPLC Column
100 X 2.1 mm, 3 µm

SLFPPPGDC20-3UM

Selectra® PFPP Guard Column
10 X 2.1 mm, 3 µm

SLGRDHLDR-HPOPT

Guard Column Holder

Introduction:

Oral fluid is a non-invasive approach to sample collection for drug screening. It is a quick procedure that is easy to collect in the field, difficult to adulterate and provides a more accurate indication of recent drug use. One of the biggest challenges in oral fluid testing is that the components of the stabilization buffer can hinder accurate data analysis and/or result in instrument downtime if not effectively cleaned up prior to injection.

This application note describes a simple and rapid solid-phase extraction (SPE) procedure for the analysis of a broad range of acidic, neutral, and basic drugs in oral fluids using UCT's Clean-Screen® XCEL I column. The functionality of the SPE sorbent effectively retains the analytes of interest while removing undesired matrix components resulting in cleaner extracts, greater method robustness and improved LC-MS/MS analysis from which both illicit and prescribed drugs can be accurately detected. HPLC separation was carried out using UCT's Selectra® PFPP column prior to detection by LC-MS/MS.

The pentafluorophenylpropyl phase can undergo dipole-dipole, and pi-pi interactions, imparting unique selectivity and retention mechanisms to the column that distinguish it from a traditional biphenyl phase. Excellent recoveries (102-129%) were obtained in all cases, except for benzodiazepines which experience losses during the wash steps at the pH employed in this procedure. It is recommended to utilize an internal standard for accurate measurement of benzophenones in this application.



Sample Pretreatment:

- To 1 mL sample (oral fluid/preservative buffer working solution), add appropriate amount of internal standard.

SPE Procedure:

1. Apply Sample Directly to SPE Column

- Load at 1 to 2 mL/minute

2. Wash Column

- 1 x 2 mL pH 6.0 Phosphate Buffer
- 1 x 3 mL DI H₂O
- 1 x 3 mL 25% Acetonitrile
- Dry column for at least 10 minutes at 80-100 psi

3. Elute (Acidic/Neutral Drugs)

- 1 x 2 mL of Hexane: Ethyl Acetate (80:20) directly into 12x75 mm culture tubes
Place tubes aside for collection of basic drugs as directed in step 5.

4. Wash Column

Note: This step will mitigate matrix effects introduced by the oral fluid.

- 1 x 2 mL of MeOH: DI H₂O (50:50)
- Dry columns for 5-10 minutes at 80-100 psi

5. Elute (Basic Drugs)

- 1 x 2 mL of Methanol: NH₄OH (98:2) directly into the same 12x75 mm culture tubes used in collection of acidic/neutral drugs.

6. Evaporate

- Evaporate the combined eluent in culture tubes to dryness
- Add 100 µL 1% HCl after 4 minutes of drying

7. Reconstitute

- Reconstitute sample in 200 µL MeOH: H₂O (50:50)



LC-MS/MS Parameters		
System	Shimadzu LC30AD w/ MS-8050	
Column	Selectra® PFPP HPLC Column 100 X 2.1 mm, 3 µm (P/N: SLPFPP100ID21-3UM)	
Guard Column	Selectra® PFPP Guard Column 10 X 2.1 mm, 3 µm (P/N: SLPFPPGDC20-3UM)	
Column Temperature	40 °C	
Column Flow Rate	0.3 mL/min	
Injection Volume	5 µL	
Autosampler Temperature	10 °C	
Gradient Program		
Time (min)	% Mobile Phase A 5 mM Ammonium Formate with 0.1% Formic Acid in H ₂ O	% Mobile Phase B 5 mM Ammonium Formate with 0.1% Formic Acid in MeOH
0	100	0
8.0	0	100
10.0	0	100
10.1	100	0
15.0	100	0




Figure 1: Chromatogram of Solvent Standard, 5 ng/mL




Figure 2: Chromatogram of Extracted Standard, 5 ng/mL






Figure 3: Hydromorphone Calibration Curve (0.5, 1, 2.5, 5, 10, 20 ng/mL) $R^2 = 0.9982$

Absolute Recovery - Oral Fluid				
Analyte	1 ng/mL (n=3)	Rel. Std Dev %	5 ng/mL (n=3)	Rel. Std Dev %
Alpha-hydroxyalprazolam	35%	5%	42%	7%
Alprazolam	129%	5%	65%	1%
Amitriptyline	120%	1%	101%	7%
Amphetamine	118%	2%	100%	5%
Atenolol	117%	6%	99%	3%
Bupivacaine	117%	5%	92%	2%
Buprenorphine	119%	3%	98%	7%
CBD	103%	3%	84%	6%
CBN	117%	2%	90%	6%
Clonazepam	118%	5%	105%	7%
Cocaine	117%	5%	91%	5%
Codeine	120%	2%	99%	7%
COOH-THC	108%	9%	86%	4%
Dextrophan	117%	2%	91%	4%
Diazepam	107%	6%	101%	6%
EDDP	117%	3%	86%	4%
EME	114%	2%	96%	2%
Fentanyl	115%	2%	94%	6%
Hydrocodone	116%	3%	94%	7%
Hydromorphone	115%	2%	98%	8%
Imipramine	119%	1%	94%	4%
Ketamine	119%	3%	94%	3%
6-MAM	113%	3%	95%	3%
Lorazepam	36%	8%	40%	7%
MDPV	116%	3%	90%	6%
Meperidine	118%	1%	93%	6%
Methamphetamine	118%	2%	94%	4%
Methylphenidate	118%	6%	94%	3%
Midazolam	118%	5%	85%	6%
MDEA	120%	5%	95%	7%
MDMA	115%	4%	93%	6%
Morphine	113%	2%	96%	5%
Naltrexone	117%	6%	98%	6%
Norbuprenorphine	109%	3%	99%	7%
Norcodeine	115%	6%	105%	9%
Nordiazepam	102%	3%	101%	2%
Norketamine	122%	2%	100%	7%
Normeperidine	112%	4%	91%	5%
Nortriptyline	117%	1%	93%	4%
OH-THC	113%	7%	90%	5%
Oxazepam	51%	7%	55%	7%
Oxycodone	117%	4%	98%	5%
Oxymorphone	113%	3%	100%	4%
PCP	104%	8%	100%	6%
Phentermine	117%	2%	99%	5%
Temazepam	47%	5%	56%	8%
THC	114%	3%	92%	9%
Tramadol	117%	4%	93%	4%
Zolpidem	114%	1%	93%	2%



Conclusion:

Use of UCT's protocol for the screening of oral fluids is efficient in the extraction of acidic, neutral, and basic drugs while minimizing matrix effects introduced by oral fluid buffer components. The use of the highly selective Clean Screen® XCEL I SPE column allows for increased analyte sensitivity, enhanced specificity for selected functional groups, low organic solvent consumption, and optimized chromatographic resolution.

In this procedure, a wash step using 3 mL DI H₂O followed by 3 mL 25% Acetonitrile was used prior to elution of analytes of interest. The choice of 25% acetonitrile was optimal for removal of dyes and buffer salts in the storage mechanism while ensuring optimal recoveries.

Also critical to this application was the use of a two-step elution procedure which allowed the wash of oral fluid buffer constituents from the column between elution of acidic/neutral, and basic drug eluents.

Excellent results were obtained for all drug classes, except for benzodiazepines which experienced loss during the wash steps at the pH employed in this procedure. Utilization of an internal standard will aid in accurate measurement of the benzodiazepine compounds studied in this application.

9110-02-01

UCT, LLC • 2731 Bartram Road • Bristol, PA 19007 800.385.3153 • 215.781.9255

www.unitedchem.com Email: methods@unitedchem.com

©UCT, LLC 2023 • All rights reserved

