


Solid Phase Extraction of Novel Synthetic 2-Benzylbenzimidazole Opioid Compounds "Nitazenes"



UCT Part Numbers

CSDAU133

Clean Screen® DAU
3 mL, 130 mg

SCS27-C18GDC21

SelectraCore® C18
Guard Column
5 x 2.1 mm, 2.7 µm

SPHPHO7001-10

Select pH buffer
pouch 100 mM
Phosphate buffer pH 7.0

SCS27-C181021

SelectraCore® C18 Column
100 x 2.1 mm, 2.7 µm

SLGRDHLDR-HPOPT

Selectra® Direct
Connect Guard Holder

Summary:

A novel group of synthetic opioids, known as benzylbenzimidazole-opioids or "nitazenes," is emerging. Originally synthesized in the 1950s as a potential analgesic, these compounds were never approved for medical use.¹ However, they are now resurfacing as a significant threat within the ongoing opioid epidemic. These potent synthetic opioids range from three to twenty times more potent than fentanyl.² The CFRSE detected the first nitazene analog, isotonitazene, in biological samples in the United States back in July 2019.³ Since then, the number of cases nationwide has continued to increase. To address this issue, the Drug Enforcement Administration (DEA) has temporarily classified eight nitazenes as schedule I substances.⁴ This application note introduces a simple targeted extraction method for the analysis of nine nitazene compounds from urine and blood utilizing UCT's flagship Clean Screen® DAU column and new SelectraCore® C18 core-shell column.



Sample Pretreatment:

- In a test tube add 0.5 mL of sample, internal standard, 200 μ L of acetonitrile (ACN), and 1.3 mL of 100 mM phosphate buffer pH 7.
- Vortex and centrifuge samples for 10 minutes at 3000 rpm.

SPE Procedure:

1. Condition Column

- a) 1 x 3 mL MeOH
- b) 1 x 3 mL phosphate buffer pH 7

2. Load Sample

- a) Load at 1 to 2 mL/minute

3. Wash Column

- a) 1 x 3 mL DI H₂O
- b) 1 x 3 mL 50:50 MeOH:H₂O

4. Dry Column

- a) Dry column for at least 10 minutes under full pressure or vacuum

5. Elute Analytes

- a) 1 x 3 mL of MeOH:NH₄OH (98:2)
- Note:** Prepare elution solvent daily

6. Evaporate

- a) Add 250 μ L of 10% HCl in methanol and vortex
- b) Evaporate eluate at 35°C, 10 psi

7. Reconstitute

- a) 1 mL of 50:50 MeOH:H₂O

Notes:

- Centrifuging blood samples causes a decrease in sample recovery, but improves visual cleanliness of sample at the end of SPE
- As an alternative to adding 10% HCl before evaporation, samples can be evaporated at 30°C, 5 psi



LC-MS/MS Parameters		
LC-MS/MS System	Shimadzu Nexera LC-30AD with MS-8050	
UHPLC Column	SelectraCore® C18 Column 100 x 2.1 mm, 2.7 μm (PN: SCS27-C181021)	
Guard Column	SelectraCore® C18 Guard Column 5 x 2.1 mm, 2.7 μm (PN: SCS27-C18GDC21)	
Column Temperature	40°C	
Flow Rate	0.45 mL/min	
Injection Volume	5 μL	
Gradient Program		
Time (min)	% Mobile Phase A: 0.1% formic acid in Water	% Mobile Phase B: 0.1% formic acid in Methanol
0	90	10
2.5-3.5	57	43
7	30	70
8-11	0	100
11.3-15	90	10

MRM Table:

Analyte	Parent Ion (m/z)	Product Ion 1 (m/z)	CE (eV)	Product Ion 2 (m/z)	CE (eV)	RT (min)
Butonitazene	425.5	100.1	-23	72.1	-45	5.83
Clonitazene	386.5	100.1	-26	125.1	-36	4.03
Etonitazene	397.4	100.1	-21	72.0	-36	3.88
Etonitazepyne	395.6	98.1	-23	56.1	-55	3.80
Flunitazene	371.3	100.1	-23	73.1	-26	3.41
Isotonitazene	411.5	100.1	-21	72.2	-45	4.53
Metodesnitazene	339.2	100.1	-21	72.1	-40	2.09
Metonitazene	383.5	100.1	-22	72.2	-39	3.38
Protonitazene	411.7	100.1	-24	72.1	-39	4.98

* CE = collision energy, RT = retention time



Chromatogram:




Figure 1: Chromatogram of extracted 15 ng/mL blood sample

Example Solvent Calibration Curves:




Figure 2a: Example of a 6-point solvent calibration curve for Clonitazene with linear equation & R^2 value (0.5, 1, 2.5, 5, 10 & 20 ng/mL)




Figure 2b: Example of a 6-point solvent calibration curve for Protonitazene with linear equation & R^2 value (0.5, 1, 2.5, 5, 10 & 20 ng/mL)

Results:

Urine (n=5)									
	1 ng/mL			5 ng/mL			15 ng/mL		
Analyte	Recovery (%)	RSD (%)	Matrix Effects (%)	Recovery (%)	RSD (%)	Matrix Effects (%)	Recovery (%)	RSD (%)	Matrix Effects (%)
Butonitazene	97	4	-14	93	2	-12	93	2	-12
Clonitazene	101	6	-5	93	2	0	95	2	-3
Etonitazene	103	6	-5	94	1	-1	99	1	-6
Etonitazepyne	104	6	-2	95	1	4	96	3	-1
Flunitazene	102	5	-5	100	4	1	98	2	-2
Isotonitazene	100	4	-4	94	2	1	98	2	-1
Metodesnitazene	98	2	8	93	4	5	106	4	-10
Metonitazene	97	3	0	92	2	3	98	1	4
Protonitazene	100	4	-5	94	0	1	96	1	-3

Blood (n=5)									
	1 ng/mL			5 ng/mL			15 ng/mL		
Analyte	Recovery (%)	RSD (%)	Matrix Effects (%)	Recovery (%)	RSD (%)	Matrix Effects (%)	Recovery (%)	RSD (%)	Matrix Effects (%)
Butonitazene	75	5	9	74	6	14	81	3	0
Clonitazene	87	3	6	83	5	12	89	3	-3
Etonitazene	93	4	0	87	5	8	94	1	-5
Etonitazepyne	94	2	13	89	4	18	97	1	0
Flunitazene	96	5	-3	92	3	5	98	2	-11
Isotonitazene	85	4	8	83	5	15	90	4	-1
Metodesnitazene	95	6	11	89	4	17	94	6	-6
Metonitazene	95	4	8	89	4	16	98	3	-1
Protonitazene	87	4	3	81	3	10	87	3	-5

Recovery was calculated by comparing peak area of pre-spiked samples to peak area of post-spiked samples. Matrix effects were calculated by comparing peak area of post-spiked samples to peak area of evaporated solvent standards. A negative matrix effect indicates ion suppression while a positive matrix effect indicates ion enhancement.



Conclusion/Discussion:

A simple extraction method was developed for the extraction of nine nitazene compounds from urine and blood. Analytes were extracted using UCT's flagship column Clean Screen® DAU and analyzed on a LC-MS/MS equipped with UCT's new SelectraCore® C18 core-shell column. All analytes were separated in 6 minutes with a short total run time of 15 minutes. Isomers protonitazene and isotonitazene were successfully separated on the core-shell column. Due to these compounds' novelty and potency, developing an extraction with a low limit of quantitation was crucial and challenging.


A sizeable amount of the non-polar analytes, particularly butonitazene and isotonitazene, remain in the test tube after loading the sample onto the SPE cartridge. To better retain the analytes, 200 µL of acetonitrile was added during sample preparation. This is vital for detection and quantitation at low concentrations. Another discovery made during method development was that free-base nitazene compounds are volatile. It was difficult to avoid the evaporation step after extraction as these compounds are present at low concentrations in biological samples. Like amphetamines, hydrochloride acid was added to the elution solvent before evaporation to create more stable salt forms.

The extraction method was evaluated using quality control samples prepared at low, medium, and high concentrations. Recovery and matrix effect for each analyte were calculated using pre-spiked samples, post-spiked samples, and evaporated solvent standards. Pre-spiked samples are extracted biological samples spiked during sample preparation. Post-spiked samples are extracted biological samples spiked after the extraction into the elution solvent. Evaporated solvent standards are spiked elution solvent samples with 10% HCl in methanol that were fully dried and reconstituted. Formulas for recovery and matrix effect are shown below:

$$\% \text{ Recovery} = \frac{\text{Peak Area of Pre-Spiked Samples}}{\text{Peak Area of Post-Spiked Samples}} \times 100$$

$$\text{Matrix Effect} = \left(\frac{\text{Peak Area of Post-Spiked Samples}}{\text{Peak Area of Evaporated Samples}} \right) \times 100$$

Extraction recoveries of analytes from urine ranged from 93-106% with relative standard deviations less than 10%. Matrix effects for all analytes were within $\pm 25\%$ making it easy to implement and validate this method in laboratories that follow ANSI/ASB Standard 036. Extraction recoveries of analytes from blood ranged from 74-96% with relative standard deviations less than 10%. Matrix effects for all analytes in blood were also within $\pm 25\%$.



References:

1. Diversion Control Division, Benzimidazole-Opioids Other Name: Nitazenes (2022).
2. Vandeputte, M.M., Krotulski, A.J., Walther, D. et al. Pharmacological evaluation and forensic case series of Npyrrolidino etonitazene (etonitazepyne), a newly emerging 2-benzylbenzimidazole 'nitazene' synthetic opioid. Arch Toxicol 96, 1845–1863 (2022).
<https://doi.org/10.1007/s00204-022-03276-4>
3. Alex J Krotulski, Donna M Papsun, Sherri L Kacinko, Barry K Logan, Isotonitazene Quantitation and Metabolite Discovery in Authentic Forensic Casework, Journal of Analytical Toxicology, Volume 44, Issue 6, July 2020, Pages 521–530, <https://doi.org/10.1093/jat/bkaa016>
4. Seven Benzimidazole-Opioids: Butonitazene, Etodesnitazene, Flunitazene, Metodesnitazene, Metonitazene, N-Pyrrolidino Etonitazene, and Protonitazene, 86 Fed. Reg. 69183-69186 (December 7, 2021)

UCT, LLC • 2731 Bartram Road • Bristol, PA 19007 • 800.385.3153 • 215.781.9255

www.unitedchem.com • Email: methods@unitedchem.com

©UCT, LLC 2023 • All rights reserved

