

# Solid Phase Extraction of Natural Cannabinoids and Metabolites from Blood and Urine Emily Eng\*, Stephanie Reichardt, Ritesh Pandya and Abderrahim Abdelkaoui | UCT, Inc.

## INTRODUCTION

The Cannabis plant contains over a hundred identified natural cannabinoids with  $\Delta^9$  - tetrahydrocannabinol ( $\Delta^9$  -THC),  $\Delta^8$  -tetrahydrocannabinol ( $\Delta^8$ -THC), and Cannabidiol (CBD) being some of the most well-known. Both marijuana and hemp are forms of *Cannabis*, only differentiated by their  $\Delta^9$ -THC content. Cannabis is considered to be marijuana when all parts of the plant, whether growing or not, have a  $\Delta^9$  -THC content above 0.3% of its dry weight [1]. It's becoming increasingly more important to both accurately identify and quantitate cannabinoids in biological matrices due to the changing legal status of marijuana across the U.S. and the rise in popularity of cannabinoids for medical and recreational use [2].


While cannabinoids and their analysis are not novel, there is much room for improvement as they continue to persist as a problem for toxicology laboratories due to their "sticky" nature. SPE procedures generally yield low recoveries and liquid-liquid extractions can be tedious for analysts to perform and generally require large amounts of solvent. This poster outlines two separate extraction methods for four natural cannabinoids and  $\Delta^9$ -THC metabolites from blood and urine using UCT's Clean Screen® THC and Styre Screen® HLB respectively. Using a Shimadzu LCMS-8050 Triple Quadrupole Mass Spectrometer, all six analytes were separated in a short convenient 12-minute method. Equipped with UCT's SelectraCore® C18 core-shell column, this LC method successfully separates isomers  $\Delta^8$  -THC and  $\Delta^9$  -THC. Both solid-phase extraction procedures can be easily implemented with good recoveries and low matrix effects into high throughput laboratories.


#### [1] 21 U.S.C. § 802 (16) (2022)

[2] Pellati, Federica et al. "Cannabis sativa L. and Nonpsychoactive Cannabinoids: Their Chemistry and Role against Oxidative Stress, Inflammation, and Cancer." BioMed research international vol.2018 1691428. 4 Dec. 2018, doi:10.1155/2018/1691428

## INSTRUMENT PARAMETERS Shimadzu Nexera LC-30AD w/ MS-8050 **LC-MS/MS System UHPLC Column** SelectraCore® C18 Column 100 x 2.1 mm, 2.7 µm **Guard Column** SelectraCore® C18 5 x 2.1 mm, 2.7 µm **Column Temperature** 40°C **Flow Rate** 0.4 mL/min **Injection Volume Mobile Phase A** 0.1% formic acid in water **Mobile Phase B** 0.1% formic acid in methanol **Gradient** Conc. B 50% (0 min) – 80% (3 min) – 90% (7.5 min) – 100% (8 to 9 min) – 50% (9.1 to 12 min)

### Chromatogram of extracted 50 ng/mL urine sample, minutes 4.10 – 7.20





# SPE PROCEDURE

Blood Extraction Clean Screen THC 6mL, 200mg (PN: CSTHC206) Urine Extraction Styre Screen HLB 3mL, 60mg (PN: SSHLB063)

• 1mL sample + ISTD + 1mL of ACN + 1ml phosphate buffer







- Vortex and centrifuge • 1 x 2mL of MeOH
- 1 x 2mL of pH 7 phosphate buffer
- 3 Load Load at 1 to 2 mL/minute
  - 1 x 3mL deionized water • 1 x 3mL 50% MeOH in deionized water
- 5 Dry • Dry column for at least 10 minutes under full pressure or vacuum
- 1 x 3mL of 60:40 MeOH: Hexane 6 Elution Note: shake or vortex elution solvent well before use
- 7 Dry Elute • Evaporate eluate under a constant gentle stream of nitrogen ≤ 40°C

y = 20075x - 61124

 $R^2 = 0.9996$ 

0 20 40 60 80 100 120


Concentration, ng/mL

• Reconstitute in 1mL of MeOH Alternative compatible solvents or volumes can be used



- 0.5mL sample + ISTD + 2mL cold ACN: Acetone (75:25) Vortex well and centrifuge
  Decant into 3mL pH 7 phosphate buffer
- 1 x 2mL of MeOH • 1 x 2mL of pH 7 phosphate buffer
- Load at 1 to 2 mL/minute
- 2 x 3mL deionized water • 2 x 3mL 40% MeOH in deionized water
- Dry column for at least 10 minutes under full pressure or vacuum
- 1 x 3mL of 89:9:2 ACN:MeOH:Acetic Acid
- Evaporate eluate under a constant gentle stream of nitrogen ≤ 40°C
- Reconstitute in 1mL of MeOH Alternative compatible solvents or volumes can be used

# CALIBRATION CURVES




2.5E+07


2.0E+07


1.5E+07


5.0E+06

₹ 1.0E+07

















Styre Screen® HLB and Clean Screen® THC were utilized to extract natural cannabinoids from urine and blood respectively. Urine extractions resulted in excellent recoveries and low matrix effects. All analytes at three different concentrations had recoveries > 90% and matrix effects within  $\pm$ 25%. Blood extractions also resulted in good recoveries with low matrix effects for most analytes. All analytes at three different concentrations had recoveries > 74% and matrix effects within  $\pm 26\%$ .

# CONCLUSION



Two different solid-phase extraction methods have been developed for the analysis of natural cannabinoids from blood and urine. A short 12-minute LC-MS/MS method using a SelectraCore $^{\circ}$  C18 core-shell column can separate isomers,  $\Delta 8$ -THC and  $\Delta 9$ -THC. Cannabinoids were extracted from urine using Styre Screen® HLB with high recoveries and low matrix effects. The addition of 1 mL of acetonitrile (ACN) during sample preparation of the urine samples was optimized to prevent analytes from sticking to the sample test tube and loads more analyte onto the SPE cartridge. This results in high recoveries. After a protein precipitation, Clean Screen® THC was used to extract analytes from blood. The extracts had good recoveries and low matrix effects for all analytes. Using ACN:acetone (75:25) as the protein precipitation solvent resulted in the best compromise of high recoveries and lower matrix effects compared to other solvents such as ACN and ACN with 0.1% formic acid.

Disclosure: The speaker, author, moderator, planning member and/or presenter/s do have financial relationships with UCT, Inc., as defined in the AACC policy on potential bias or conflict of interest. The specific product/s: Styre Screen® HLB, Clean Screen® THC and the SelectraCore® C18 column will be mentioned and/or discussed.





UCT, Inc. 2731 Bartram Road | Bristol, PA 19007 Phone: (800)385-3153 Email: info@unitedchem.com URL: unitedchem.com



