



Analysis of Pesticide and Mycotoxin Residues in Cannabis using QuEChERS Extraction, ChloroFiltr® dSPE Cleanup and LC-MS/MS

Abderrahim Abdelkaoui | UCT, Inc. 2731 Bartram Road, Bristol, PA 19007



INTRODUCTION

With the recent trends in legalization, interest in cannabis and cannabis-based products (e.g. concentrated oils, soda, candy and other edible forms) have dramatically increased. Pesticide and mycotoxin residues can pose significant health risks, especially with chronic exposure. The warm, wet conditions that are ideal for growing cannabis are also conducive to the growth of molds and fungi which can produce carcinogenic mycotoxins, including ochratoxin A and aflatoxins. As a result, testing for the presence of pesticides and mycotoxins in cannabis is essential to ensure consumer safety. With the widespread legalization of cannabis, UCT is presenting this simple method which would be beneficial for any research or production facility wanting to implement regulatory testing.

LC-MS/MS PARAMETERS




HPLC Column	UCT Selectra® PFPP, 100 × 2.1 mm, 3 µm (p/n: SLFPFP100ID21-3UM)
Guard column	UCT Selectra® PFPP, 10 × 2.0 mm, 3 µm (p/n: SLFPFPGDC20-3UM)
Guard column holder	UCT Selectra® Guard Cartridge Holder (p/n: SLGRDHLDR)
Column temperature	40°C
Flow rate	0.400 µL/min
Injection volume	2 µL

GRADIENT PROGRAM

TIME (min)	Mobile Phase A (%) 10 mM Ammonium Formate with 0.1% Formic Acid in DI Water	Mobile Phase B (%) Acetonitrile
0.0	98	2
12.0	0	100
13.0	0	100
13.1	98	2
16.5	98	2

MS CONDITIONS

MS/MS system	ABSCIEX QTrap 6500+
Ionization Mode	Electrospray Ionization in positive mode (ESI+)
Ion Spray Voltage (IS)	+4500.00
Temperature (TEM)	300°C
Curtain Gas (CUR)	40
Ion Source Gas 1 (GS1)	50
Ion Source Gas 2 (GS2)	50



Calibration Curve Examples - Spirotetramat, Dimethomorph E and Ochratoxin A
7 point calibration curve prepared at 0.5, 1, 5, 10, 25, 50 and 100 ng/mL.



UCT, Inc.
2731 Bartram Road | Bristol, PA 19007
Phone: (800)385-3153
Email: info@unitedchem.com



unitedchem.com



Disclosure: The speaker, author, moderator, planning member and/or presenter/s do have financial relationships with UCT, Inc., as defined in the AACCP policy on potential bias or conflict of interest. The specific product/s: Selectra® PFPP column, SpinFiltr®, ChloroFiltr® and UCT QuEChERS will be mentioned and/or discussed.

© 2023 UCT, Inc. All Rights Reserved. | Questions: methods@unitedchem.com




Figure 1: Critical Separation of Dimethomorph Isomers




Figure 2: Separation of Critical Isomers - Spinosyn A and D




Figure 3: Retention of Early Eluting Analytes - Daminozide and Acephate




Figure 4: Pesticides Panel TIC

Extract Cannabis Flower

Sample Preparation:

- 10 grams of cannabis flower was ground and homogenized using 2 mL DI water.
- The sample was mixed well in a Spex 2010 Geno/Grinder® for 10 min.
- The sample was thoroughly mixed and vortexed to achieve homogeneity.
- 10 different samples were weighed at 1 gram each.
- 5 samples were spiked at low (5 ng) and 5 samples at high (25 ng) fortification levels.

Extraction Procedure:

- Place each prepared sample in a 50 mL centrifuge tube.
- Add 10 µL of Internal Standard(s).
- Add 5 mL of DI water to each sample and vortex mix well to ensure the analyte concentration is distributed as equally as possible throughout sample.
- Add 10 mL of acetonitrile containing 2% formic acid.
- Add the contents of the ECMSSC-MP Mylar pouch (4 g MgSO₄ and 1 g NaCl) and shake for 10 minutes using the Spex 2010 Geno/Grinder®.
- The sample is centrifuged at ≥ 3000 × g for 5 minutes.

Cleanup Procedure:

- Transfer 1 mL aliquot of supernatant into ECQSF154CT dSPE cleanup tube containing 50 mg MgSO₄, 150 mg Endcapped C18, 150 mg ChloroFiltr® and 150 mg PSA.
- Vortex the sample for 30 seconds.
- Centrifuge the sample at ≥ 3000 × g for 5 minutes.
- Transfer the purified and filtered sample extract into an autosampler vial for analysis on ABSciex 6500+ Triple Quad LC-MS/MS.

CONCLUSION

This Poster outlines a QuEChERS method for the simultaneous analysis of cannabis for 67 pesticides and 5 mycotoxins residues in cannabis flower. Sample purification is carried out using UCT's new cleanup product SpinFiltr®, which combines the convenience of classical dispersive-SPE (dSPE) with an ultrafiltration tube containing a 0.2 µm filter membrane to simultaneously remove unwanted matrix components and filter the sample prior to LC analysis.

The SpinFiltr® dSPE tube uses MgSO₄, PSA, C18 and ChloroFiltr® for sample cleanup. ChloroFiltr® is a unique sorbent designed for the removal of chlorophyll and unlike graphitized carbon black (GCB), does not result in the loss of planar analytes. Liquid chromatography, using a Selectra® PFPP column, coupled to tandem mass spectrometry (LC-MS/MS) is used for analysis of the pesticides. The injection time is 17 minutes with 11 minutes of scan time. The method achieved good separation of all analytes in the panel, especially some critical isomers as Dimethomorph E and Z, Spinosyn A and D, and Mevinphos E and Z.

The method was evaluated by fortifying cannabis samples with each compound at low and high concentrations (n=5 each). The average recovery obtained was predominantly in the range of 70-100% and the RSD ≤20%. For some analytes, lower recoveries were obtained due to polarity, as in the case of Daminozide or sensitivity, as in the case of Cyfluthrin.