

# Simultaneous Analysis of 19 Novel Synthetic Cannabinoids in Urine Using SPE and LC-MS/MS



## UCT Part Numbers

**SSHLD063**

Styre Screen® HLD  
60 mg, 3mL Column

**SLC-18100ID21-18UM**

Selectra® C18 HPLC Column  
100 X 2.1 mm, 1.8 µm

**SPPHO6001-10**

Select PH Buffer Pouches 100 mM  
Phosphate pH 6.0

**SLGRDHLDR-HPOPT**

UHPLC Direct Connect  
Guard Holder

**SLC-18GDC20-18UMOPT**

Selectra® C18 Guard Column  
10 X 2.0 mm, 1.8 µm

## Summary:

Newly identified synthetic cannabinoids pose a significant threat to public health and safety, as their implications in drug overdose and adverse events continue to rise in the United States and around the world. The diverse chemical structures of these compounds have a great impact on their potency and side effects. These synthetic cannabinoids were previously un-reported in forensic toxicology casework in the United States. There are currently few published methods available for the analysis of these novel compounds. However, the importance of identifying and extracting these compounds from various biological matrices is becoming more critical for accurate forensic criminal investigations and clinical diagnostics.

This application note outlines a solid-phase extraction (SPE) and LC-MS/MS method for the analysis of 19 synthetic cannabinoids in urine. These specific compounds were selected based on positivity rates from several key testing labs in the area. The use of UCT's Styre Screen® HLD highly crosslinked polymeric SPE sorbent ensures efficient extraction of the synthetic cannabinoids while removing undesired matrix components and yielding clean extracts. LC separation was carried out using a Selectra® C18 UHPLC column which resulted in excellent retention and baseline separation of the critical isobaric compounds ADBICA N-pentanoic acid and ADB-PINACA N-pentanoic acid metabolite in under 10 minutes.



## SPE Procedure:

### 1. Sample Preparation

- To 1 mL of urine add 1 mL of pH 6 phosphate buffer (0.1M) and internal standard(s)
- Mix/vortex briefly

**Note:** A hydrolysis protocol may be required if conjugated compounds are to be included into the above drug panel

### 2. Condition Cartridge

- 1 x 1 mL MeOH
- 1 x 1 mL DI H<sub>2</sub>O

### 3. Apply Sample

- Load sample at 1- 2 mL/minute

### 4. Wash Cartridge

- 1 x 2 mL DI H<sub>2</sub>O
- 1 x 2 mL ACN: H<sub>2</sub>O (20:80, v/v) containing 1% Formic Acid
- Dry cartridges under full vacuum or pressure for 5 minutes

### 5. Elute Analytes

- 1 x 3 mL Ethyl Acetate
- Collect at 1-2 mL/minute

### 6. Dry Eluate

- Evaporate to dryness at < 40°C.


### 7. Reconstitute

- Reconstitute sample in 1 mL of mobile phase or other appropriate organic solvent.



| LC-MS/MS Parameters |  |   |
|---------------------|--|---|
| LC-MS/MS System     | Shimadzu LCMS-8050   |   |
| UHPLC Column        | Selectra® C18 (100 X 2.1 mm, 1.8 µm)                         |   |
| Guard Column        | Selectra® C18 (10 X 2.0 mm, 1.8 µm)                          |   |
| Column Temperature  | 50°C   |   |
| Flow Rate           | 0.4 mL/min   |   |
| Injection volume    | 5 µL   |   |
| Gradient Program    |  |   |
| Time (min)          | % Mobile Phase A:<br>0.1% formic acid in DI H <sub>2</sub> O | % Mobile Phase B:<br>0.1% formic acid in MeOH |
| 0                   | 70   | 30  |
| 3.5                 | 55   | 35  |
| 9                   | 0  | 100   |
| 11                  | 0  | 100   |
| 11.1                | 70   | 30  |
| 14                  | 70   | 30  |

## Chromatogram:



**Figure 1:** Chromatogram of 25 ng/mL extracted sample demonstrating the isobaric separation of ADBICA N-pentanoic acid and ADB-PINACA N-pentanoic acid metabolite.




## Results:

| Recovery (% , n=5)                         |           |         |          |         |
|--|-----------|---------|----------|---------|
| Analyte                                    | 2.5 ng/mL | RSD (%) | 10 ng/mL | RSD (%) |
| AB-PINACA N-Pentanoic Acid Metabolite      | 102%      | 0.16    | 110%     | 0.32    |
| ADBICA N-Pentanoic Acid                    | 115%      | 0.12    | 108%     | 0.71    |
| ADB-PINACA N-Pentanoic Acid Metabolite     | 87%       | 0.12    | 113%     | 1.04    |
| AB-FUBINACA Oxobutanoic Acid               | 92%       | 0.06    | 97%      | 0.14    |
| 5-Fluoro ADBICA                            | 91%       | 0.08    | 106%     | 0.17    |
| ADB-BICA                                   | 88%       | 0.09    | 101%     | 0.19    |
| 4-cyano CUMYL-BUTINACA                     | 105%      | 0.09    | 112%     | 0.28    |
| ADB-FUBICA                                 | 86%       | 0.12    | 106%     | 0.13    |
| 5-Fluoro MDMB-PICA                         | 100%      | 0.11    | 109%     | 0.25    |
| PB-22 3-Carboxyindole Metabolite           | 97%       | 0.12    | 107%     | 0.18    |
| MDMB-FUBICA                                | 98%       | 0.09    | 108%     | 0.19    |
| BB-22-Carboxyindole Metabolite             | 103%      | 0.06    | 106%     | 0.18    |
| UR-144 (XLR11) N-Pentanoic Acid Metabolite | 95%       | 0.11    | 105%     | 0.20    |
| AKB-48 N-Pentanoic Acid Metabolite         | 100%      | 0.09    | 109%     | 0.12    |
| MDMB-FUBICA                                | 93%       | 0.24    | 95%      | 0.30    |
| AB-CHMINACA 3-methyl Butanoic Acid         | 99%       | 0.11    | 110%     | 0.24    |
| BB-22                                      | 93%       | 0.27    | 103%     | 0.66    |
| MA-CHMINACA                                | 95%       | 0.38    | 97%      | 0.49    |
| MDMB-CHMINACA                              | 100%      | 0.43    | 100%     | 1.02    |



## Representative Calibration Curve (ADBICA N-Pentanoic Acid):



## Conclusions:

This application note outlines a simple SPE procedure for the analysis of 19 synthetic cannabinoids in urine using UCT's Styre Screen® HLD highly crosslinked polymeric SPE cartridges. All 19 compounds were analyzed in under 10 minutes using LC-MS/MS. The use of a Selectra® C18 UHPLC column resulted in excellent peak shape for all the compounds included in the method, including baseline separation of any isobaric compounds. The recoveries obtained in this research study are satisfactory for the vast majority of analytes despite their diverse chemical structures. Extracted urine samples fortified at two concentrations (2.5 and 10 ng/mL) had, on average, recoveries in the range of 85-110% and corresponding RSD values less than 5%. The quality control concentrations of 2.5 ng/mL and 10 ng/mL were chosen to ensure low-level accurate detection based on heightened potency for this class of compounds at exceptionally low biological levels. This method will be beneficial to any lab looking to implement testing of these novel synthetic cannabinoids.



**0210-01-01**

UCT, LLC • 2731 Bartram Road • Bristol, PA 19007 800.385.3153 • 215.781.9255

[www.unitedchem.com](http://www.unitedchem.com) Email: [methods@unitedchem.com](mailto:methods@unitedchem.com)

©UCT, LLC 2023 • All rights reserved

