



Determination of Diquat and Paraquat in Drinking Water by Solid-Phase Extraction and LC-MS/MS Detection



Diquat



Paraquat

UCT Part Numbers

Enviro-Clean® RFV0050CT
50 mL centrifuge tubes


Enviro-Clean® SPE cartridge:
EUCCX11Z
Carboxylic acid 100 mg/10 mL

Summary:

Diquat and paraquat are fast-acting, non-selective herbicides used widely as desiccants and defoliants. They are quaternary amines that are highly water soluble. Their toxicity and presence in bodies of water adversely affect aquatic life and human health. Therefore, it is essential to determine their levels in drinking water samples.

The traditional drinking water method for diquat and paraquat analysis is EPA method 549.2. This method employs an ion-pairing reverse phase (C8) Solid-Phase Extraction (SPE) followed by ion-pairing HPLC with UV or photodiode array detection. The traditional method is time-consuming (extracting 250 mL sample), needs ion-pairing reagents, and is less sensitive than alternative extraction and analysis options.

This application outlines a novel weak cation exchange SPE method with LC-MS/MS detection for diquat and paraquat. The method is fast and sensitive, using only 10 mL of water sample. In addition, there is no need for ion-pairing. Moreover, quaternary amines are retained onto the sorbent by a cation exchange mechanism; washing the sorbent with organic solvents after extraction will not wash off the retained amines; however, it will provide a much cleaner extract than using the traditional reverse phase C8 sorbent.



Notes

Diquat and paraquat cations absorb onto glass surfaces. Plastic labware was used for the entire procedure.

Deuterated diquat and paraquat are not stable in aqueous solutions and were added to the final extracts as instrumental internal standards.

Preparation of buffers, elution solvent, and mobile phase:

A. 400 mM phosphate buffer (pH 7)

Dissolve 20.9 g of potassium phosphate dibasic and 10.9 g of potassium phosphate monobasic in 500 mL of reagent water. Adjust pH to 7 with diluted potassium hydroxide or phosphoric acid.

B. 25 mM phosphate buffer (pH 7)

Mix 50 mL of solution A. with 750 mL reagent water.

C. 25 mM ammonium formate buffer (pH 8)

Weigh 1.6 g of ammonium formate to a 1-L volumetric flask, add 950 mL reagent water and 1.4 mL of ammonium hydroxide and mix well. Adjust pH to 8 with diluted formic acid or ammonium hydroxide. Dilute to mark with reagent water.

D. Elution solvent: 10% formic acid in acetonitrile

Mix 10 mL of formic acid with 90 mL of acetonitrile (MeCN), and mix well.

E. Mobile phase buffer: 100 mM ammonium acetate buffer (pH 5)


Weigh 7.78 g of ammonium acetate and 2 g of glacial acetic acid into a 1-L mobile phase reservoir, and add 998 mL of reagent water. Sonicate for 30 min to dissolve the salt and acid, and remove the dissolved gases.

Sample pretreatment:

Transfer 10 mL of water sample to a 50 mL centrifuge tube (RFV0050CT), add 25 μ L of 400 mM phosphate buffer (pH7), and spike with appropriate amounts of diquat and paraquat standards for fortified samples, cap and mix well.

SPE Procedure:


- 1) Place the labeled SPE cartridges (EUCCX11Z) onto the glass block manifold lid.
- 2) Condition the cartridges with 3 mL of methanol (MeOH), and 3 mL of 25 mM phosphate buffer (pH 7).
- 3) Load the pretreated water samples onto the SPE cartridges, and apply a low vacuum for a slow dropwise flow (about 2-3 mL/min).
- 4) Wash the 50 mL centrifuge tubes with 3 mL of 25 mM ammonium formate buffer (pH 8) and apply the rinsate to the cartridges. Repeat with 3 mL of MeOH.
- 5) Dry the cartridges by applying full vacuum for 3 min.
- 6) Insert labeled 12*75 mm polypropylene test tubes into the manifold.
- 7) Elute with 3*1 mL of 10% formic acid in MeCN, pass 1/3 through, soak for 1 min, and draw the remaining through slowly.
- 8) Evaporate the eluates to dryness under a stream of nitrogen in a 45°C water bath.
- 9) Reconstitute with 900 μ L of the mobile phase (100 mM ammonium acetate buffer (pH5): MeCN, 30:70, v/v), add 100 μ L of 1 ppm IS mix, vortex and transfer 200 μ L to 250- μ L polypropylene inserts held in 2-mL vials.
- 10) Extracts are ready for analysis.



LC-MS/MS Method	
HPLC: Thermo Scientific Dionex UltiMate 3000® LC System	
Column: Thermo Scientific, Acclaim® Trinity™ Q1, 50 x 2.1 mm, 3 µm	
Guard Column: Thermo Scientific, Acclaim® Trinity™ Q1, 10 x 2.1 mm, 3 µm	
Column Temperature: 25 °C	
Column Flow Rate: 0.300 mL/min	
Auto-sampler Temperature: 10 °C	
Injection Volume: 5 µL	
Mobile phase (isocratic): 30% of 100 mM ammonium acetate buffer (pH 5) and 70% of MeCN	

MS Parameters	
Polarity	ESI +
Spray voltage V	3500 V
Vaporizer Temperature	400 °C
Ion transfer capillary temperature	350 °C
Sheath gas pressure	30 arbitrary units
Auxiliary gas pressure	15 arbitrary units
Q1 and Q3 peak width (FWHM)	0.4 and 0.7 Da
Collision gas and pressure	Ar at 2.3 mTorr
Scan type	SRM
Cycle time	1 sec
Acquisition method	EZ Method

SRM transitions							
Compound	Rt (min)	Precursor ion	Product ion 1	CE 1	Product ion 2	CE 2	S-lens (V)
Paraquat d8	2.31	192.94	177.60	24	164.71	30	53
Paraquat	2.33	184.95	168.52	17	114.66	23	59
Diquat d3	3.52	185.97	157.65	22	130.61	34	55
Diquat	3.53	183.08	156.63	22	129.63	33	55




Results:



Recovery and RSD% Obtained from 6 Replicated Fortified Water Samples

Compound	Spiked at 0.5 µg/L		Spiked at 25 µg/L	
	Recovery %	RSD % (n=6)	Recovery %	RSD % (n=6)
Paraquat	96.1	7.1	97.9	5.2
Diquat	89.2	7.0	87.9	7.1

Chromatogram of a Water Sample Fortified with 0.5 µg/L of Diquat and Paraquat



Matrix Matched Calibration Curves (Dynamic Linearity Range: 0.5 – 50 µg/L)



DCN-312280-283

UCT, LLC • 2731 Bartram Road • Bristol, PA 19007 800.385.3153 • 215.781.9255

www.unitedchem.com Email: methods@unitedchem.com

©UCT, LLC 2013 • All rights reserved

