


LSD and Metabolites In Blood, Plasma/Serum, and Urine Using: 200 mg Clean Screen® Extraction Column



UCT Part Numbers

ZSDAU020
Clean Screen® DAU
200 mg, 10 mL Column

Procedure:

1. Prepare Sample

- To 1 mL of 100 mM phosphate buffer (pH 6) add internal standards.*
- Add 1 mL of whole blood, serum/ plasma, urine.
- Add 2 mL of 100 mM phosphate buffer (pH 6).
- Vortex and centrifuge as appropriate.

2. Condition Column

- 1 x 3 mL CH₃OH.
 - 1 x 3 mL D.I. H₂O.
 - 1 x 1 mL 100 mM phosphate buffer (pH 6).
- Note:** aspirate at < 3 inches Hg to prevent sorbent drying out

3. Apply Sample

- Load sample at 1-2 mL / minute.

4. Wash Column

- 1 x 3 mL DI H₂O.
- 1 x 3 mL 100 mM acetic acid.
- 1 x 3 mL CH₃OH.
- Dry column (5 minutes at > 10 inches Hg).

5. Elute LSD and Metabolites

- 1 x 3 mL ethyl acetate; acetonitrile: ammonia (78: 20: 2 v/v)
- OR
- 1 x 3 mL CH₂Cl₂/ IPA/ ammonia (78:20:2 v/v)
 - Collect eluate at 1-2 mL /minute.

6. Evaporation

- Evaporate eluates to dryness under a gentle stream of nitrogen.

7. Reconstitute

- Reconstitute sample in 50 µL of CH₃OH.
- Inject 5 µL.




Instrument Conditions		
Column	100 x 2.1 mm (3 µm) Selectra® Phenyl (UCT, LLC)	
Flowrate	0.5 mL/ minute Column	
Temperature	Ambient Detector: API 3200 QTRAP MS/MS	
Mobile Phase		
Time (min)	Acetonitrile %	0.1% Formic acid
0	30	70
3	90	10
3.1	30	70
5.0	30	70

Compound	MRM Transistion	Cerilliant #
LSD	324.2/ 223.1	-005
Iso-LSD	324.2/281 (223.1)	I-010
Nor-LSD	310.2/209.1	L-017
OH-LSD	356.2/338.1	O-013
*LSD-D3	327.2/226.1	L-002



Chromatogram Of LSD And Metabolites



Recovery: > 90% (N=10)

LOD: 0.1 ng/mL

DCN-903020-77

UCT, LLC • 2731 Bartram Road • Bristol, PA 19007 800.385.3153 • 215.781.9255

www.unitedchem.com Email: methods@unitedchem.com

©UCT, LLC 2023 • All rights reserved

