




FORENSICS



UCT

sample prep  
made easy



GLASS BLOCK MANIFOLD

## 16-Port and 24-Port Vacuum Manifold

Description	Units	Part Number	Part Number
Complete Vacuum Manifold	1	16 Position System VMF016GL	24 Position System VMF024GL

NOTE: A Complete Vacuum Manifold System includes a glass block, Corian® manifold lid, a cover gasket, a vacuum gauge and valve assembly, teflon tips (16 or 24), an adjustable collection rack, bulkhead luer fittings, plugs (16 or 24), and manifold safety tray.


## Vacuum Manifold Parts Accessories

Description	Units	Part Number	Part Number
		16 Position System	24 Position System
Glass Block	1	VMF06123	VMF04123
Manifold Lid	1	VMF06120	VMF04120
Gasket	2	VMF06121	VMF04121
Collection Rack	1	VMF06125	VMF04125
Vacuum Gauge & Valve Assembly	1	VMF02122	
Teflon Tips	12	VMF020TT	
Retaining Clips for Collection Rack	12	VMF02129	
Bulkhead Luer Fittings	12	VMF21BFN	
Manifold Lid Legs	4	VMF02120-1	
Luer Lock Plugs	12	VMF21PLN	
Clean-Thru® Tips	50	CLTTP050	
Manifold Safety Tray	1	VMF02072	
Adaptors for 1mL, 3mL, 6mL, 10mL or 15mL column	15	AD0000AS	
Stopcocks:	16	VMF02116	
	24	VMF02024	
Vacuum Pump	1	PMPV192A	


- 1/8 hp, 115 volts, 4.2 amps, 60 Hz



Glass Block



16 Position Lid



24 Position Lid



Gasket



Collection Racks



Vacuum Gauge & Valve Assembly



Teflon Tips



Retaining Clips for Collection Rack



Bulkhead Luer fittings



Manifold Lid Legs



Luer-Lock Plugs



Clean-Thru Tips



Manifold Safety Tray



Stopcocks



Adaptors



Vacuum Pump

## PRICES AND TERMS

Our prices are subject to change without notice. The price in effect when we receive your order will apply. All prices are in US Dollars and are F.O.B. Lewistown, PA 17044. Terms of payment are net 30 days.

## MINIMUM ORDERS

We welcome all orders, therefore, we do not have a minimum order requirement. When ordering, please include your purchase order number, complete "Ship To" and "Bill To" address, catalog number, quantity, and description of product(s). Also include your name and a phone number where you can be reached should we have any questions concerning your order.

## SHIPMENTS

Normal processing is within 24 hours after receipt of an order. Unless special shipping requests have been made, our trained staff will send all orders by UPS Ground service. The appropriate shipping charges (freight & insurance costs) will be added to the invoice, unless otherwise instructed by the customer.

## SPECIAL PRICING

We offer special pricing for volume purchases and standing orders. These discounts apply to bonded phase extraction column purchases only. Please call a sales representative for more information on special pricing qualifications.

## RETURN POLICY

Our Quality Manager will handle all returns. Before returning merchandise, please call to obtain a return authorization number from the Quality Manager. We will need to know the reason for the return, date of purchase, purchase order number and invoice number in order to issue a return authorization number. Return merchandise must be received before a credit can be issued. Returns will not be accepted after 90 days. A restocking fee of 25% of the price paid, or a minimum of \$25.00 (whichever is greater) will be charged on all returns.

## WARRANTY

All products manufactured by UCT are guaranteed against defects in materials and workmanship for a period of 90 days after shipment. UCT will replace any items that prove to be defective during this time period.

The exclusive remedy requires the end user to first advise UCT of the defective product by phone or in writing. Secondly, the defective product must be returned within 30 days after proper approval from our Quality Manager. All returns must indicate the purchase order number, the lot number and the shipping date. UCT's total liability is limited to the replacement cost of UCT products.

This warranty does not apply to damage resulting from misuse.

### Placing An Order

Phone: 215.781.9255  
800.541.0559  
Fax: 215.785.1226

### Mail

UCT, LLC  
2731 Bartram Rd.  
Bristol, PA 19007

### Technical Support

Phone: 215.781.9255  
800.385.3153  
Fax: 215.785.1226  
Email: [info@unitedchem.com](mailto:info@unitedchem.com)  
Forum: [forums.unitedchem.com](http://forums.unitedchem.com)  
Web: [www.unitedchem.com](http://www.unitedchem.com)

**25** CELEBRATING  
**YEARS**  
1986 \* 2011

DCN-112270-13



FORENSICS