




UNITED CHEMICAL TECHNOLOGIES, INC.

Technical Data Sheet

Product name: A0300
Chemical name: ACETOXYETHYLTRICHLOROSILANE
Formula: C₄H₇Cl₃O₂
Mol. Weight: 221.5
Cas No.: 0018204-80-3

Structure:



Selling specifications

<u>Tests</u>	<u>Specifications</u>	<u>Test methods</u>
IR	MUST MATCH REFERENCE	UCT0050
PURITY (SUM OF MAIN PEAK + ONE ISOMER)	97 % Min.	UCT0376
COLOR	Report only (GARDNER)	UCT0032
CUSTOMER SHELF LIFE:	6	Months from the date of shipment

*The specification values are only valid in conjunction with the methods listed.
The shelf life is only valid if the material is stored in the original unopened container.*

2731 Bartram Road / Bristol PA 19007 USA / Phone 215: 781-9255 / Fax 215: 785-1226
Email: info@unitedchem.com Internet: [HTTP://www.unitedchem.com](http://www.unitedchem.com)