



UNITED CHEMICAL TECHNOLOGIES, INC.

Technical Data Sheet

Product name: A0394.5
Chemical name: ACRYLOXYPROPYLEMETHYLDIMETHOXYSILANE
Formula: C₉H₁₈O₄Si
Mol. Weight: 236.1
Cas No.: 0013732-00-8

Selling specifications

Tests

IR
PURITY
CUSTOMER SHELF LIFE:

Specifications

MUST MATCH REFERENCE
93.0 % Min.
6

Test methods

UCT0050
UCT0219
Months from the date of shipment

The specification values are only valid in conjunction with the methods listed.

The shelf life is only valid if the material is stored in the original unopened container.

2731 Bartram Road / Bristol PA 19007 USA / Phone 215: 781-9255 / Fax 215: 785-1226
Email: info@unitedchem.com Internet: [HTTP://www.unitedchem.com](http://www.unitedchem.com)