




UNITED CHEMICAL TECHNOLOGIES, INC.

Technical Data Sheet

Product name: B2500
Chemical name: BIS(TRIMETHYLSILYL)ACETAMIDE
Common name: BSA
Mol. Weight: 203.40

Description: LIQUID

Structure:



Selling specifications

Tests

PURITY

IR

CUSTOMER SHELF LIFE:

Specifications

92 % Min.

MUST MATCH REFERENCE

6

Test methods

UCT0333

UCT0050

Months from the date of shipment

Physical and chemical properties

SP. GRAVITY @20DEG.C

BOILING POINT

REF. INDEX N(20,D)

FLASH POINT

0.83

71 Deg.C.

1.42

17 Deg.C.

ASTM D891-95

DIN 5175

UCT0044

UCT0381

*The specification values are only valid in conjunction with the methods listed. These methods are available upon request.
The shelf life is only valid if the material is stored in the original unopened container.*