




UNITED CHEMICAL TECHNOLOGIES, INC.

Technical Data Sheet

Product name: D6010
Chemical name: DIPHENYLDIMETHOXYSILANE
Formula: C₁₄H₁₆O₂Si
Mol. Weight: 244.40
Cas No.: 0006843-66-9

Description: CLEAR, COLORLESS LIQUID

Structure:



Selling specifications

Tests

PURITY
IR
CUSTOMER SHELF LIFE: 6

Specifications

97 % Min.
MUST MATCH REFERENCE

Test methods

UCT0368
UCT0050
Months from the date of shipment

Physical and chemical properties

SP. GRAVITY @20DEG.C
REF. INDEX N(20,D)
FLASH POINT
BOILING POINT (@15 MM)

1.08
1.5447
250 Deg.C.
161 Deg.C.

ASTM D891-95
UCT0044
UCT0381
DIN 5175

*The specification values are only valid in conjunction with the methods listed.
The shelf life is only valid if the material is stored in the original unopened container.*