



# UNITED CHEMICAL TECHNOLOGIES, INC.

## *Technical Data Sheet*

**Product name:** PS072

**Chemical name:**

DIMETHYLSILOXANE ETHYLENE OXIDE PROPYLENE OXIDE COPOLYMER (60% propylene oxide, 40% ethylene oxide)

**Cas No.:** 0067762-85-0

## *Selling specifications*

### Tests

IR  
VISCOSITY @ 25 DEG.C  
SP. GRAVITY @25DEG.C  
REF. INDEX N(25,D)  
**CUSTOMER SHELF LIFE:** 12

### Specifications

MUST MATCH REFERENCE  
1500 to 2000 cS  
1.0165 to 1.0372  
1.4435 to 1.4485

### Test methods

UCT0050  
UCT0031  
ASTM D891-95  
UCT0044  
Months from the date of shipment

*The specification values are only valid in conjunction with the methods listed.*

*The shelf life is only valid if the material is stored in the original unopened container.*

2731 Bartram Road / Bristol PA 19007 USA / Phone 215: 781-9255 / Fax 215: 785-1226  
Email: [info@unitedchem.com](mailto:info@unitedchem.com) Internet: [HTTP://www.unitedchem.com](http://www.unitedchem.com)