



UNITED CHEMICAL TECHNOLOGIES, INC.

Technical Data Sheet

Product name: PS073
Chemical name: DIMETHYLSILOXANE-ETHYLENE OXIDE COPOLYMER
Cas No.: 0068937-57-2

Selling specifications

Tests

SP. GRAVITY @20DEG.C
REF. INDEX N(25,D)
COLOR
VISCOSITY @ 25 DEG.C
CUSTOMER SHELF LIFE: 12

Specifications

Report only
Report only
Report only (GARDNER)
80 to 150 cS

Test methods

ASTM D891-95
UCT0044
UCT0032
UCT0035
Months from the date of shipment

The specification values are only valid in conjunction with the methods listed.

The shelf life is only valid if the material is stored in the original unopened container.

2731 Bartram Road / Bristol PA 19007 USA / Phone 215: 781-9255 / Fax 215: 785-1226
Email: info@unitedchem.com Internet: [HTTP://www.unitedchem.com](http://www.unitedchem.com)