



UNITED CHEMICAL TECHNOLOGIES, INC.

Technical Data Sheet

Product name: PS200
Chemical name: GLASSCLAD 18

Description: Dark liquid

Selling specifications

Tests

STABILITY TEST 10% WATER SOLN. MUST BE CLEAR AFTER 30 MIN.
APPEARANCE FREE OF PARTICULATES
COLOR 11 GARDNER Max.
CUSTOMER SHELF LIFE: 6

Specifications

Test methods

UCT0062
UCT0063
UCT0032
Months from the date of manufacturing

Physical and chemical properties

FLASH POINT 10 Deg.C.
SP. GRAVITY @25DEG.C 0.88
SOLIDIFICATION POINT -30 Deg.C.

UCT0381
ASTM D891-95
PROPERTY

PS200 is a proprietary monomeric octadecylsilane derivative

The specification values are only valid in conjunction with the methods listed.

The shelf life is only valid if the material is stored in the original unopened container.

2731 Bartram Road / Bristol PA 19007 USA / Phone 215: 781-9255 / Fax 215: 785-1226
Email: info@unitedchem.com Internet: [HTTP://www.unitedchem.com](http://www.unitedchem.com)