



# UCT SPECIALTIES, LLC

## *Technical Data Sheet*

**Product name:** PS225  
**Chemical name:** GLASSCLAD EG  
**Common name:** GLASSCLAD EG  
**Cas No.:** NA

**Description:** Liquid

## *Selling specifications*

### Tests

IR  
SOLIDS  
SP. GRAVITY @25DEG.C  
**CUSTOMER SHELF LIFE:**

### Specifications

MUST MATCH REFERENCE  
18 to 22 %  
Report only  
3

### Test methods

UCT0050  
UCT0033  
ASTM D891-95  
Months from the date of shipment

Also called GLASSCLAD EG

*The specification values are only valid in conjunction with the methods listed.  
The shelf life is only valid if the material is stored in the original unopened container.*

2731 Bartram Road / Bristol PA 19007 USA / Phone 215: 781-9255 / Fax 215: 785-1226  
Email: [info@unitedchem.com](mailto:info@unitedchem.com) Internet: [HTTP://www.unitedchem.com](http://www.unitedchem.com)