



UNITED CHEMICAL TECHNOLOGIES, INC.

Technical Data Sheet

Product name: PS252
Chemical name: GLASSCLAD FF
Cas No.: NA

Selling specifications

Tests

VISCOSITY @ 25 DEG.C
COLOR
SOLIDS
SP. GRAVITY @25DEG.C
SKIN OVER
TACK FREE
IR
CUSTOMER SHELF LIFE: 2

Specifications

350 to 600 cS
10 APHA Max.
48.0 to 52.0 %
.9200 to 1.000
Report only (MIN)
Report only (MIN)
MUST MATCH REFERENCE

Test methods

UCT0031
UCT0032
UCT0033
ASTM D891-95
UCT0068
UCT0068
UCT0050
Months from the date of shipment

*The specification values are only valid in conjunction with the methods listed.
The shelf life is only valid if the material is stored in the original unopened container.*

2731 Bartram Road / Bristol PA 19007 USA / Phone 215: 781-9255 / Fax 215: 785-1226
Email: info@unitedchem.com Internet: [HTTP://www.unitedchem.com](http://www.unitedchem.com)