



UNITED CHEMICAL TECHNOLOGIES, INC.

Technical Data Sheet

Product name: T2845
Chemical name: (3,3,3-TRIFLUOROPROPYL)TRICHLOROSILANE
Formula: C₃H₄Cl₃F₃Si
Mol. Weight: 231.50
Cas No.: NA

Selling specifications

Tests

IR
PURITY

CUSTOMER SHELF LIFE:

Specifications

MUST MATCH REFERENCE
97 % Min.

6

Test methods

UCT0050
UCT0488

Months from the date of shipment

*The specification values are only valid in conjunction with the methods listed.
The shelf life is only valid if the material is stored in the original unopened container.*

2731 Bartram Road / Bristol PA 19007 USA / Phone 215: 781-9255 / Fax 215: 785-1226
Email: info@unitedchem.com Internet: [HTTP://www.unitedchem.com](http://www.unitedchem.com)