




UNITED CHEMICAL TECHNOLOGIES, INC.

Technical Data Sheet

Product name: T2951
Chemical name: TRIMETHYLCHLOROSILANE (99%)
Formula: C₃H₉ClSi
Mol. Weight: 108.70

Description: LIQUID

Structure:



Selling specifications

Tests

PURITY
IR
COLOR
CUSTOMER SHELF LIFE: 6

Specifications

99 % Min.
MUST MATCH REFERENCE
50 APHA Max.

Test methods

UCT0013
UCT0050
UCT0032
Months from the date of shipment

Physical and chemical properties

SP. GRAVITY @20DEG.C
BOILING POINT
REF. INDEX N(20,D)
FLASH POINT

0.86
57 Deg.C.
1.3890
-18 Deg.C.

ASTM D891-95
DIN 5175
UCT0044
UCT0381

*The specification values are only valid in conjunction with the methods listed.
The shelf life is only valid if the material is stored in the original unopened container.*