



# UCT Specialties, LLC

## Technical Data Sheet

**Product Name:** A0701-KG

**Chemical:** N (2- AMINOETHYL) 3-AMINOPROPYLTRIMETHOXYSILANE TECHNICAL GRADE

**Cas Number:** 1760-24-3

Test	Specifications	Method
Appearance Description	Colorless to straw colored liquid	SOP82-4.080
IR	Pass	SOP82-4.075
Base Neutral Equivalent by Titration	110.0000 - 120.0000	SOP82-4.148
Refractive Index N(25,D)	1.4420 - 1.4480	SOP82-4.070

*The specification values are only valid in conjunction with the methods listed.*

*The shelf life for this product is 6 months from the date of shipment.*

*The shelf life is only valid if the material is stored in the original unopened container.*

*The specifications and lot analysis are only valid using the listed methods.*

*The test methods and specifications were current when the product was lotted.*