



# UCT Specialties, LLC

## Technical Data Sheet

**Product Name:** A0800-KG  
**Chemical:** 3-AMINOPROPYLTRIMETHOXY SILANE  
**Cas Number:** 13822-56-5

Test	Specifications	Method
Purity by GC	97.0000 Min.	SOP82-4.062
IR	Pass	SOP82-4.075
IR Search Results	0.9000 Min.	SOP82-4.075

*The specification values are only valid in conjunction with the methods listed.*  
*The shelf life for this product is 6 months from the date of shipment.*  
*The shelf life is only valid if the material is stored in the original unopened container.*  
*The specifications and lot analysis are only valid using the listed methods.*  
*The test methods and specifications were current when the product was lotted.*