



UCT Specialties, LLC

Technical Data Sheet

Product Name: PC088.3-KG
Chemical: PLATINUM OCTANAL/OCTANOL COMPLEX
Cas Number: 68412-56-6

Test	Specifications	Method
Appearance Description	Amber solution, no particulates	SOP82-4.080
Cure Test	25.0000 Min.	SOP82-4.145
ICP of Pt	2.0000 - 2.2000	SOP82-4.277

*The specification values are only valid in conjunction with the methods listed.
The shelf life for this product is 6 months from the date of shipment.
The shelf life is only valid if the material is stored in the original unopened container.
The specifications and lot analysis are only valid using the listed methods.
The test methods and specifications were current when the product was lotted.*