



UCT Specialties, LLC

Technical Data Sheet

Product Name: PS401-KG

Chemical: POLYDIMETHYLSILOXANE, AMINOALKYL T-STRUCTURE (80-130 cS)

Cas Number: 71750-80-6

Test	Specifications	Method
Appearance	Pass	VISUAL
IR	Pass	SOP82-4.075
Solids @ 150 Deg. C		SOP82-4.065
Viscosity in cS	80.0000 - 130.0000	SOP82-4.063
Refractive Index N(25,D)	1.4025 - 1.4125	SOP82-4.070
Specific Gravity @ 25 Deg. C	0.9650 - 0.9850	SOP82-4.207
Amine Value	0.2500 - 0.7500	SOP82-4.078

The specification values are only valid in conjunction with the methods listed.

The shelf life for this product is 12 months from the date of shipment.

The shelf life is only valid if the material is stored in the original unopened container.

The specifications and lot analysis are only valid using the listed methods.

The test methods and specifications were current when the product was lotted.