



UCT Specialties, LLC

Technical Data Sheet

Product Name: PS426-KG

Chemical: 1.0%-VINYLMETHYL, DIMETHYLSILOXANE COPOLYMER, TRIMETHYL
TERMINATED (1000 cS)

Cas Number: 67762-94-1

Test	Specifications	Method
IR	Pass	SOP82-4.075
Methylvinylsiloxane (Wt.%)	0.7500 Min.	SOP82-4.081
Specific Gravity @ 25 Deg. C	Report Only	SOP82-4.207
Viscosity in cS	900.0000 - 1100.0000	SOP82-4.063

The specification values are only valid in conjunction with the methods listed.

The shelf life for this product is 12 months from the date of shipment.

The shelf life is only valid if the material is stored in the original unopened container.

The specifications and lot analysis are only valid using the listed methods.

The test methods and specifications were current when the product was lotted.