



ENVIRO



UCT

ChloroFiltr®

UCT
ChloroFiltr
ECMPSSGG15CT
MgSO₄/PSA/CHF

UCT
ChloroFiltr
ECMPSSGG15CT
MgSO₄/PSA/CHF

UCT
ChloroFiltr
ECMPSSGG15CT

Selectively removing
Chlorophyll without
Compromising your
Planar Pesticides

INNOVATION THROUGH CHEMISTRY






ChloroFiltr®



ChloroFiltr® is a polymeric sorbent available exclusively from UCT. It is designed to replace graphitized carbon black (GCB) for the efficient removal of chlorophyll without loss of planar analytes.

ChloroFiltr® has been tested against hundreds of pesticides and herbicides and has been shown to reduce chlorophyll concentration by greater than 82% without loss of planar analytes.



Don't Sacrifice Your Planar Pesticides!



6106-01-01

800.385.3153

www.unitedchem.com