



# UCT Specialties, LLC

## Technical Data Sheet

**Product Name:** D5200-KG

**Chemical:** (N,N-DIMETHYL-3-AMINOPROPYL)TRIMETHOXYSILANE

**Cas Number:** 2530-86-1

<b>Test</b>	<b>Specifications</b>	<b>Method</b>
Appearance	Pass	VISUAL
IR	Pass	SOP82-4.075
Purity by GC	95.0000 Min.	SOP82-4.062

*The specification values are only valid in conjunction with the methods listed.*

*The shelf life for this product is 6 months from the date of shipment.*

*The shelf life is only valid if the material is stored in the original unopened container.*

*The specifications and lot analysis are only valid using the listed methods.*

*The test methods and specifications were current when the product was lotted.*