



UCT Specialties, LLC

Technical Data Sheet

Product Name: D5636-KG
Chemical: DIMETHYL-n-OCTADECYLCHLOROSILANE
Cas Number: 18643-08-8

Test	Specifications	Method
Purity (AS ALKYLTRICHLOROSILANES)		SOP82-4.108
Purity (AS ALKYLTRICHLOROSILANES)	97.0000 Min.	SOP82-4.108
IR	Pass	SOP82-4.075
Color APHA		SOP82-4.064

*The specification values are only valid in conjunction with the methods listed.
The shelf life for this product is 6 months from the date of shipment.
The shelf life is only valid if the material is stored in the original unopened container.
The specifications and lot analysis are only valid using the listed methods.
The test methods and specifications were current when the product was lotted.*