



# UCT Specialties, LLC

## Technical Data Sheet

**Product Name:** E6250-KG

**Chemical:** 2-(3,4-EPOXYCYCLOHEXYL)ETHYLTRIMETHOXYLANE

**Cas Number:** 3388-04-3

Test	Specifications	Method
Appearance	Pass	VISUAL
Purity (AS SUM OF ALPHA AND BETA ISOMERS)	97.0000 Min.	SOP82-4.108
IR	Pass	SOP82-4.075
IR Search Results	0.9000 Min.	SOP82-4.075

*The specification values are only valid in conjunction with the methods listed.*

*The shelf life for this product is 6 months from the date of shipment.*

*The shelf life is only valid if the material is stored in the original unopened container.*

*The specifications and lot analysis are only valid using the listed methods.*

*The test methods and specifications were current when the product was lotted.*