



# UCT Specialties, LLC

## Technical Data Sheet

**Product Name:** H6860-KG

**Chemical:** (HEPTADEC AFLUORO-1,1,2,2-TETRAHYDRODECYL) TRICHLOROSILANE

**Cas Number:** 78560-44-8

Test	Specifications	Method
Purity by GC	97.0000 Min.	SOP82-4.062
IR	Pass	SOP82-4.075
Appearance	Pass	VISUAL
IR Search Results	0.9000 Min.	SOP82-4.075

*The specification values are only valid in conjunction with the methods listed.*

*The shelf life for this product is 6 months from the date of shipment.*

*The shelf life is only valid if the material is stored in the original unopened container.*

*The specifications and lot analysis are only valid using the listed methods.*

*The test methods and specifications were current when the product was lotted.*