



# UCT Specialties, LLC

## Technical Data Sheet

**Product Name:** M9160-KG

**Chemical:** N-METHYL-N-TRIMETHYLSILYLTRIFLUOROACETAMIDE

**Cas Number:** 24589-78-4

Test	Specifications	Method
Appearance	Pass	SOP82-4.080
Purity by GC	97.0000 Min.	SOP82-4.062
IR	Pass	SOP82-4.075

*The specification values are only valid in conjunction with the methods listed.*

*The shelf life for this product is 24 months from the date of shipment.*

*The shelf life is only valid if the material is stored in the original unopened container.*

*The specifications and lot analysis are only valid using the listed methods.*

*The test methods and specifications were current when the product was lotted.*