



UCT Specialties, LLC

Technical Data Sheet

Product Name: O9780-KG
Chemical: n-OCTADECYLTRIMETHOXY SILANE
Cas Number: 3069-42-9

Test	Specifications	Method
H+ NMR		SOP82-4.239
Chloride		SOP82-4.069
Purity (SUM OF ISOMERS) by GC		SOP82-4.108
Appearance	Pass	SOP82-4.080
Purity (SUM OF ISOMERS) by GC	96.0000 Min.	SOP82-4.108
IR	Pass	SOP82-4.075

*The specification values are only valid in conjunction with the methods listed.
The shelf life for this product is 6 months from the date of shipment.
The shelf life is only valid if the material is stored in the original unopened container.
The specifications and lot analysis are only valid using the listed methods.
The test methods and specifications were current when the product was lotted.*