



UCT Specialties, LLC

Technical Data Sheet

Product Name: O9818.3-KG
Chemical: 7-OCT-1-ENYLTRIMETHOXY-SILANE
Cas Number: 52217-57-9

Test	Specifications	Method
Appearance	Pass	VISUAL
C13 NMR	Pass	SOP82-4.239
Proton NMR - R&D	Pass	SOP82-4.239
Si NMR	Pass	SOP82-4.239
Purity (SUM OF ISOMERS) by GC	97.0000 Min.	SOP82-4.108
IR	Pass	SOP82-4.075
IR Search Results	0.9000 Min.	SOP82-4.075

*The specification values are only valid in conjunction with the methods listed.
The shelf life for this product is 6 months from the date of shipment.
The shelf life is only valid if the material is stored in the original unopened container.
The specifications and lot analysis are only valid using the listed methods.
The test methods and specifications were current when the product was lotted.*