



UCT Specialties, LLC

Technical Data Sheet

Product Name: PS034-KG

Chemical: POLYDIMETHYLSILOXANE, TRIMETHYLSILOXY TERMINATED

Cas Number: 63148-62-9

Test	Specifications	Method
Appearance	Pass	VISUAL
Purity by GC	97.0000 Min.	SOP82-4.062
IR	Pass	SOP82-4.075
Refractive Index N(25,D)	1.3700 - 1.3800	SOP82-4.070
Specific Gravity @ 25 Deg. C	0.7510 - 0.7710	SOP82-4.207
Viscosity in cS	0.5800 - 0.7800	SOP82-4.299
Color APHA		SOP82-4.064

The specification values are only valid in conjunction with the methods listed.

The shelf life for this product is 12 months from the date of shipment.

The shelf life is only valid if the material is stored in the original unopened container.

The specifications and lot analysis are only valid using the listed methods.

The test methods and specifications were current when the product was lotted.