



# UCT Specialties, LLC

## Technical Data Sheet

**Product Name:** PS076-KG

**Chemical:** (N-TRIMETHOXYSILYLPROPYL)POLYETHYLENIMINE, 50 % in IPA (135-200 cS)

**Cas Number:** 136856-91-2

<b>Test</b>	<b>Specifications</b>	<b>Method</b>
Appearance	Pass	VISUAL
IR	Pass	SOP82-4.075
Viscosity in cS	135.0000 - 200.0000	SOP82-4.063
Solids @ 100 Deg. C	46.5000 - 53.5000	SOP82-4.065
Amine Value	4.5000 - 6.3000	SOP82-4.296
Specific Gravity @ 25 Deg. C	0.8900 - 0.9300	SOP82-4.207
Refractive Index N(25,D)	1.4300 - 1.4400	SOP82-4.070

*The specification values are only valid in conjunction with the methods listed.*

*The shelf life for this product is 12 months from the date of shipment.*

*The shelf life is only valid if the material is stored in the original unopened container.*

*The specifications and lot analysis are only valid using the listed methods.*

*The test methods and specifications were current when the product was lotted.*