



# UCT Specialties, LLC

## Technical Data Sheet

**Product Name:** S1588-KG  
**Chemical:** STYRYLETHYLTRIMETHOXY-SILANE  
**Cas Number:** 134000-44-5

Test	Specifications	Method
Purity (SUM OF ISOMERS) by GC		SOP82-4.108
Appearance	Pass	SOP82-4.080
H+ NMR	Report Only	SOP82-4.239
IR	Pass	SOP82-4.075
Purity (ALKYL & ALKENYLPHENETHYLTRIMETHOXY-SILANES)	97.0000 Min.	SOP82-4.108
IR Search Results	0.9000 Min.	SOP82-4.075

*The specification values are only valid in conjunction with the methods listed.  
The shelf life for this product is 6 months from the date of shipment.  
The shelf life is only valid if the material is stored in the original unopened container.  
The specifications and lot analysis are only valid using the listed methods.  
The test methods and specifications were current when the product was lotted.*